

Vetlux Vet Products Global Supplier of Veterinary Consumables





"Special acknowledgment to our parent company Vetmed.kz for all the support and assistance."

#### ABOUT US

Founded in March 2024, Vetlux embodies the spirit of renewal and exploration much like the season itself, when swallows return and nature awakens to new possibilities. As a subsidiary of Vetmed, we are deeply grateful for the foundational support, product expertise, and market insights provided by our parent company. Without Vetmed, there would be no Vetfux.

Under the visionary leadership of Max, head of Vetmed, Vetlux benefits from his exceptional expertise and dedication. As a highly experienced entrepreneur and frands-on leader, Max brings profound knowledge in veterinary products and a genuine passion for pet medical devices and consumables. His guidance has been instrumental in shaping our vision and operations, and we are immensely thankful for his unwavering support and mentorship.

Directly managed by Jack Wong—an adventurous and challenge-driven professional with six years of experience in the veterinary consumables market—Vetlux is driven by passion and expertise. Our VP, Kenny, is a renowned veterinary surgeon whose extensive clinical experience has earned him the trust of countless pet owners and greatly contributed to our product development and market adaptation.

Every product in our catalog undergoes rigorous testing in partner clinics, and most are also validated by Vetmed before market release. With a background rooted in a human medical hospital with over 20 years of history, Vetlux possesses strong supplier relationships and refined product selection capabilities. We are the only company in China dedicated exclusively to supplying veterinary consumables.

In just one year, Vetlux has secured exclusive distribution rights in 10 countries worldwide. We strictly maintain a single-distributor policy per country to protect the interests and investments of our partners.

Beyond consumables, we also supply veterinary equipment approved by both our partner clinics and Vetmed.

Committed to social responsibility, Vetlux donates 1% of every order value in the client's name to national animal protection organizations.

At the heart of Vetlux & Vetrned lies a profound love for animals. Despite operating with minimal profit margins—sometimes even at cost—we remain unwavering in our mission. Every step forward, whether big or small, represents progress. We are dedicated to becoming a comprehensive solution for all surgical consumables in veterinary clinics, striving to make procedures safer and animals more cared-for and loved.

Though young, Vetlux possesses boundless potential and approaches every task with dedication and responsibility. We sincerely thank every customer who supports our journey. We are dreamers, and together, we are committed to making this dream grow even greater.

Driven by expertise, integrity, and a global vision, Vetlux is your trusted partner in advancing animal health and care.

"Committed to social responsibility, Vetlux donates 1% of every order value in the client's name to national animal protection organizations."

——— Max, Jack, Kenny





# Suture

"Precise, Smooth, Trusted,"

### Suture:

Original imported thread from a renowned Korean brand.

## **Needle Options:**

- Mid-range domestically produced (China) needle
- 2. Premium European-made needle
- 3. High-end MANI (Japan) needle





Structure: Braided Color: Violet or Undryed Feautre: Absorbable Code: PGA

Material:Polydioxanone



Structure: Monofilament Color: Violet Feautre: Absorbable Code: PDO

Material: Chromic Catgut



Structure: Twisted Color: Brown Feautre: Absorbable Code: CC

Material: Nylon



Structure: Monofilament Color: Black or Blue Feautre: Non-absorbable Code: NL





Structure: Twisted Color: Yellow Feautre: Absorbable Code: PC

Material:Polypropylene



Structure: Monofilament Color: Blue Feautre: Non-absorbable Code: PM

#### Material: Polvester



Structure: Braided Color: Green or White Feautre: Non-absorbable Code: PB

#### Material:Silk



Structure: Braided Color: Black or Blue Feautre: Non-absorbable Code: SK

#### Absorbable

(PGA)Polyglycolic acid (PGAR)Rapide Polyglycolic acid (PGLA)Polyglactin 910 (PDO or PDS)Polydioxanone (CC)Chromic catgut

#### Non-absorbable

(SK)Silk braided

(PB)Polyester braided

(NL)Nylon monofilament

(PM)Polypropylene monofilament

(SW)Stainless steel wire

#### Thread Length:

75cm or 90cm only for Vet

#### Thread Diameter USP:

11/0, 10/0, 9/0, 8/0, 7/0, 6/0, 5/0, 4/0, 3/0, 2/0, 0, 1, 2, 3, 4, 5

#### Needle Lenght(mm):

12mm-30mm


#### Needle Curvature:


Straight, 1/2 circle, 1/2 circle (double), 1/4 circle, 1/4 circle (double), 3/8 circle, 3/8 circle (double), 5/8 circle, Loop Round

#### Cross Section:

Round Bodied, Curved Cutting, Reverse Cutting













# Pet Skin stapler

"Fast. Precise. Healing."





Easy to use for fast, safe, and precise wound closing.

Causes minimal skin damage and helps wounds heal better.

Leaves neat stitches with barely visible scars.

Reduces procedure time and patient discomfort.

Saves time for medical staff.

Item No.	Version	Staple diamter	Nail molding size	Nail molding height
VL-38	35W	0.6	10.2	3.5
VL-39	35R	0.6	12.1	2.4





# Vet Tissue Glue

"Fast. Strong. Versatile."

### Formula:

n-butyl cyanoacrylate (similar to 3M™ Vetbond™)

#### Use:

Rapid wound closure in veterinary procedures.

#### Benefits:

Replaces sutures/bandages; bonds in seconds upon contact with moisture.

### Applications:

Minor wounds, cat claws, suture sealing, extractions, and oral surgery.



VET SKIN ADHESIVE			
Content:	3.0ml/PCS		
Features :	Drop mouth applicator for well control.		

Topical Skin Adhesive

\*Sterile



VET SKIN ADHESIVE			
Content :	0.5ml/1.0ml PCS		
Features :	Designed dropper for well coating.		













1.Remove the cap.

2.Attach the applicator and push it into place.

3.Point the tip of the applicator at the patient and apply a small amount of tissue adhesive to the surface in a thin layer. Applying too thickly will cause the adhesive to crack prematurely, bringing the skin edges together for 5-10 seconds. Applicator tip avoid contact with wet surfaces and do not apply pet biological tissue glue inside or under the epidermis.

4.After the procedure is complete, check for bleeding. If necessary, re-drop the adhesive on the edge of the wound.

#### Application:

Ideal for a variety of procedures, including: animals sterilizing operation, cat declaws, dew claw removal, closing minor lacerations, tooth extractions, and oral surgery. Also used for the adhesion of small windows, abrasions, cuts and other non-chronic wounds.









Residual Toe Excision

Lacerations

General Surgery

Eye Surgery

Sterilization Surgery

- 1. Eye Surgery: Repair of superficial damage to the eyelid or corneal.
- 2. Otorology: ear trauma repair surgery, etc.
- 3. Orthopedics: residual toe resection, fracture repair, etc.
- 4. Urology: Pet bladder stone surgery, etc.

- 5. Oncology Department: Skin turnor resection.
- 6. General Surgery: sterilization, castration, pet breast tumor resection, laceration abrasion, hernia repair, incision closure after internal fluxation of fractures, etc.
- 7. Emergency Surgery: Traumatic skin laceration.
- 8. Animal experiments: Scenarios such as back incisions of experimental animals such as mice and rats, tumor tissue block fixation, etc.





## Cat Catheter

"Flexible, kink-resistant, smooth"



#### 1.Medical-Grade Material

Lubrizol TPU construction (biocompatible) with hydrophilic coating

## 2.Low-Trauma Design

Rounded end hole + ultra-smooth surface minimize tissue damage

## 3.Enhanced Functionality

Balanced hardness for stone removal + optimized elasticity

### 4.Secure Placement

Integrated Luer connectors allow skin suture fixation





## **Canine Dual Catheter**

"Precise. Smooth. Trusted."



## Packing:

1 pc/bag, 10 bags/box; Ethylene oxide sterilized.

### Product Features:

- 1. Made of imported highly elastic medical-grade silicone for the balloon, reducing burst risk.
- 2. Features a skin-friendly soft silicone tube body and tip for enhanced comfort.
- Equipped with a fine double-lumen catheter suitable for small animals, avoiding secondary trauma from suturing and providing better care for pets.








nti backflow valve

# **Urine Collection Bag**

"Leak-proof. anti-reflux. universal-fit."

Material	Color	Size
PVC	Clear	150/300ml



Safe, compact, and comfortable for small animals.

## Secure & Adjustable:

Elastic strap ensures a stable fit for various sizes.

### Anti-Reflux Valve:

Prevents urine backflow (keep valve flat for optimal function).

## Universal Luer Connector:

Compatible with Vetlux catheters (for cats and dogs) and extension tubes.

## Easy Drainage:

Includes extension tube for hygienic urine collection.





# Self-Adhesive Bandage

"Tough. Tasteless. Trusted."

## Available in 3 Types:

1. Original:

Standard cohesive bandage.

2. Printed:

Fun patterned designs for a cheerful look.

3. Bitter:

Infused with a safe, bitter agent to deter licking and chewing.











## **Cohesive Bandage**

"Product Details "









## Easy Hand-Tear Design:

Conveniently tears by hand without the need for scissors.

## 100% Self-Adhesive:

Adheres securely to itself without requiring clips, pins, or tape.

## Fur-Friendly:

Does not stick to animal hair or fur, ensuring pain-free removal.

## Secure & Supportive:

Provides reliable, non-slip, lightweight support as the final layer of a bandage.















# **Endotracheal Tube**

"Secure, Breathe, Revive,"

## Medical-grade PVC.

For emergency airway access and tracheotomy.

Smooth tip, secure placement,

post-intubation comfort.

Standard 15mm connector.

## Medical-grade PVC

Spring-supported to prevent kinking and collapse
Bite-resistant and flexible for secure placement
Smooth cuff, standard 15mm connector, easy to intubate

## Standard

Cuffed Endotracheal Tube

ID20mm/ID25mm/ID3.0mm/ID35mm/ID4.0mm/ID4.5mm/ ID5.0mm/ID55mm/ID6.0mm/ID65mm/ID75mm/ ID8.0mm/ID85mm/ID9.0mm/ID9.5mm/ID10.0mm

Cuffed Endotracheal Tube

I.D.10.5mm/I.D.11.0mm/I.D.12.0mm/I.D.13.0mm






### Reinforced

Cuffed Endotracheal Tube

ID.2.0mm/ID.2.5mm/ID.3.0mm/ID.3.5mm/ID.4.0mm/ID.4.5mm/ID.5.0mm/ID.5.5mm/ID.6.0mm/ID.6.5mm/ID.7.5mm/ID.8.0mm/ID.8.5mm/ID.9.0mm/ID.9.5mm/ID.9.0mm/ID.9.5mm/ID.9.0mm/ID.9.5mm/ID.







# **Bird Endotracheal Tube**

"Secure, Breathe, Revive, Sizes: 1.0-6.0 Material: Soft medical silicone Features: 1. Smooth tip, radiopaque marker 2. Flexible & biocompatible Use: Anesthesia & intubation in birds



# Steel ET Tube Rack

"Sturdy. Versatile. Hygienic."









## **Vet Suction Tube**

"Sterile. Versatile. Ready."

### Features:

Medical-grade transparent PVC, non-toxic

## Full size range:

F6-F24

EO sterilized, individually packed Disposable, easy to use

## Applicable for:

Difficult sputum removal
Postoperative mucus blockage
Swallowing disorders
Oral thick secretions





## **Vet Suction Bottle**

"Sterile, Versatile, Clear,"

#### Sizes:

Fr6 to Fr18

## Packaging:

1pc/PE bag, 100pcs/carton







## Features:

Medical-grade PVC, non-toxic and transparent

## Complete set:

bottle, catheter, tube, and connector Available in 25ml/40ml volumes, some with filter





# Anesthesia Mask Set

"Secure. Gentle. Fit."

For: Dogs, cats, small animals

Uses: Anesthesia, ventilation, aerosol therapy



Transparency PC Cone Stainless Ring

PP Ring

\_





PP Ring

\_



Item No.	Size	Hight	ID	OD	Connector
VL-40	1#	44mm	24mm	48mm	15mm
VL-41	2#	49mm	28mm	53mm	IOIIIII
VL-42	3#	58mm	35mm	64mm	
VL-43	4#	98mm	29mm	92mm	22mm
VL-44	5#	132mm	49mm	123mm	2211111
VL-45	6#	141mm	71mm	136mm	

## Features:

- 1. High-transparency PC cone
- 2. Soft silicone seal
- 3. 6 sizes (#1-#6)
- 4. 15mm/22mm standard connectors







# **Manual Resuscitator**

"Versatile, Safe, Reliable,"



### Material:

Medical PVC (single-use) / Medical Silicone (reusable) Sizes: Small, Medium, Large, X-Large Animals

## **Bag Capacity:**

- 1. Resuscitation Bag: 250ml / 500ml / 1000ml / 1600ml
- 2. Reservoir Bag: Optional 600-2500ml

Connections: 15mm/22mm valve + 2.1m oxygen tube

Safety: Pressure relief valve (40/60cm H<sub>2</sub>O)



# **3way Stopcock**

"Flexible, Precise, Secure,"



Essential for fluid management in veterinary procedures.



Available in standard, 10cm, or 30cm extension models.



## Features:

- 1. Precise flow control
- 2. Leak-proof design
- 3. Flexible extension options
- 4. Easy to operate





# **Aerosol Chamber**

"Gentle, Effective, Essential,"

#### **INSTRUCTONS**

Direction for Use







Helps deliver medication for pet respiratory conditions.

#### Features:

- 1. Antistatic PP chamber
- 2. Soft silicone mask (3 sizes)

3. Safe breathing valve





1. Remove the back piece mask and the adapter. Do not remove valves during cleaning.



2. Soak both parts in lukewarm water with liquid deterg-ent for 15 minutes. Agitate gently



3. Rinse in clean















6. Replace back piece when unit is completely



# Laryngoscope Set

"Precise. Bright. Durable."

- 1. Right-curved blade design for easier and smoother intubation
- 2. Textured anti-slip handle available in Medium and Large sizes
- 3. Bright LED light transmitted through flexible fiberoptic strand
- 4. Made of medical-grade stainless steel with robust welding
- 5. Blade is fully autoclavable for repeated sterilization







# Soda Lime

"Safe, Efficient, Reliable,"

## Prevents Hypercapnia:

Safely removes CO2

## Color Change:

White to violet when saturated

## Size:

5L / bottle



## Applications:

Vet clinics, hospitals, aquariums

## **Dual Use:**

For veterinary and medical anesthesia systems





# **Vet APNEA Alarm**

"Alert. Monitor. Protect."

Easy connection to tube & machine

Alarms after 15s breathing stop

Auto-restarts if breathing resumes







Touch activation, auto power-off
USB rechargeable, sound/light alerts
Surface clean only (no autoclave)





## Anesthesia Moisture Trap

"Safe. Compatible. Reliable."



## ETCO<sub>2</sub> Water Trap

- 1. Compatible with Mindray anesthesia & ICU equipment
- 2. For respiratory accessories
- 3. Made of medical plastic (SPU)
- 4. Moisture protection with sample port



DRYLINE GAS SAMPLING LINE



# **Vet HME Filter**

"Filtering. Moisturizing. Protecting."

Filters >99.9% of bacteria/viruses



Retains heat and moisture (HME)

Fits most vet and human machines (22mm/15mm)

For anesthesia and ventilators





# **Gas Filter Canister**

"Safe, Efficient, Reliable,"

### Available Sizes:

250g: For low-flow anesthesia or small clinics

800g: For high-flow systems or large hospitals

## Replacement Weight:

250g → Replace at ≥300g

800g → Replace at ≥900g



### Maintenance:

Check every 2-6 months (varies by size/usage)

Replace if weight increases 10-20% or odor is detected





## Vet Silicone Test Lung

"Test. Calibrate. Protect."

#### Function:

Tests and calibrates veterinary anesthesia machines for safe and accurate ventilation.

Capacity	1L (adjustable for small/large animals)
Use	Machine testing, leak checks, training
Compatibility	Standard 22mm veterinary circuits
Materia	Medical-grade silicone (wipe-clean)

## Material:

Durable silicone

#### Size:

27×10cm (approx.)

### Connector:

Standard 15mm

# Package: Bag + Box

### Note:

Slight size/color variation may occur.





## Anesthesia Breathing Bag

"Soft, Efficient, Safe,"



## Disposable Latex Bags Blue, PE/PVC port, 3M tape.

Sizes: 0.5/1/2/3L.



## Reusable Silicone Bags

Off-white.

Sizes: 0.5/1/2/3L.



## Disposable Neoprene Bags

Green, neoprene or economy neoprene with PE/PVC port, 3M tape.

Sizes: 0.5/1/2/3L.



## Reusable Latex Bags

Black.

Sizes: 0.25 to 30L (10 options).





## **Breathing Coaxial Circuit**

"Flexible, Secure, Sterile,"

Length:

1.6/1.8M

Material:

Medical Grade EVA

Type:

Large Animal (22mm ID) / Small Animal (15mm ID)

Features:

Lightweight, flexible, easy to clean and sterilize

Usage:


For anesthesia machines and oxygen delivery during surgery





# Y-Type Breathing Circuit

oxygen support during surgery





## **Double Swivel Catheter Mount**

"Flexible, Secure, Sterile,"



Retainable flip-top cap (22M/15F)

7.6mm probe port + Luer gas sampling ports

Dual tabs: for bronchoscope & suction access

Maintains seal during procedures








## Vet Premium IV Cannula

"Soft. Sharp. Safe."



Size & Color: 18G–26G (color-coded for easy identification)






## Veterinary IV Catheter Needle Size

Gauge	Catheter Length	Outer diameter	Flow Rate	Color Code	Shelf Box / Shipping Case
26G	15	0.6	12	Purple •	
24G	19	0.7	22	Yellow •	
22G	25	0.9	36	Blue •	
20G	32	1	58	Pink •	100/2000
18G	32	1.2	90	Green •	
17G	45	1.5	120	White	
16G	45	1.8	172	Grey	
14G	45	2	208	Orange 🔴	



## **HOU TO USE**



Prior to venipuncture

After the needle has been advanced into the vein

The catheter has been advanced into the vein





# **Vet IV Cannula**

"Soft. Sharp. Safe."





## Sterilization:

100% EO gas sterilized



#### Sizes:

14G-24G (color-coded: green/blue/pink/white)

### Material:

Precision PET/PU catheter for kink-resistant flow

### Needle:

Ultra-sharp bevel for smooth insertion

## Packaging:

Individually sealed for sterility and convenience

### Use:

Designed for small animals and pets; single-use only





# IV Catheter - Pen type

"Soft. Sharp. Safe."



Ultra-sharp needle for minimal pain and bleeding

EO gas sterilized, individually packed

Safe and single-use for animal patients

Protective cap and vacuum sealed

Smooth catheter placement





# **Vet Extension Tube**

"Flexible. Secure. Sterile."







## Pet Infusion Set

"Flexible. Precise. Ready."

### Material:

Soft, kink-resistant medical-grade PVC

## Length:

1.5m or 3m

#### Features:

Without needle + Luer connector
(direct to catheter or disposable needle)
Ready-to-use, flexible, and durable
Precision flow regulator (standard/horned shape)



## IV Set:

Plastic needle + 50-drop filter chamber + 45PE regulator + 280cm tube + Y-port + 20cm tube + Luer lock & cap



# Infusion Set - 150ML Burette

"Sterile, Precise, Versatile,"

## Vet 150ML Burette Infusion Set

Air filter
Soft tube
Drip chamber
Flow regulator
Fiuld filter
Clamp
Stop valve
Y injection site
Protective cap







## Infusion Set - T-Connector

"Versatile, Secure, Sterile,"

#### Multi-Port Options:

1/2/3/4 Way configurations available

#### EO Sterilized:

Ethylene oxide sterilized for safety

#### Precision Filtration:

Optional 0.2µm, 1.2µm, or 5µm filters

#### Safe & Sterile:

Non-toxic, non-pyrogenic, single-use. DEHP-Free/PVC-Free tubing

#### Leak-Proof:

100% airtight and fluid-tight, blocks microorganisms

### Secure Connection:

Male/female Luer lock ends, leak-free junctions







# **Pet Infusion Pump**

"Reliable, Versatile, Precise,"

#### **PRODUCT FEATURES**

- Sound alarm: automatic infusion, sound alarm, AC/DC dual purpose.
- Quick release function: Quick release function, power on self check.
- Multiple infusion methods: infusion methods: milliliters per hour, drops per minute, and time mode.
- · Flow accuracy calibration: forget to alarm function, flow accuracy calibration.
- · Accumulated display and reset: With functions of accumulated display and reset.
- Medicine library function: memory function, medicine library function, switch between Chinese and English.
- Display and volume can be set: real-time time display and settings, and button tone settings can be adjusted.







Infusion precision



Audible & visual



Flow Control

#### PRODUCT SPECIFICATION

A.3-Inch monochrome LCD screen		
1mL/h - 1200mL/h , 0.1 m L/h		
Continuous operation using a		
new fully Charged battery		
discharge time of 12 hours		
±5%(general infusion device);Drosophila error:20%		
0.1mL - 9999.9mL		
Minimum step: 0.1mL		
Speed mode, drop count mode, time mode		
11.1V 2000 mA hAC lithium battery		
Approx. 151 * 181.5 * 205mm (L x W x H)		
About 1.5 kg		







# Pet Injection Pump

"Reliable. Intelligent. Versatile."



#### PRODUCT FEATURES

- · Adjust the alarm volume to three levels: high, medium, and low.
- · Sensitive alarm reminder ensures real-time injection safety.
- A comprehensive drug library provides medical staff with convenient operation and ensures the accuracy of injection drugs.
- The data storage and upload function facilitates medical staff to review historical alarm information.
- . Can automatically identify syringe specifications of 10 ml, 20 ml, 30 ml, and 50 ml.



#### PRODUCT SPECIFICATION


Syringe size	10mL, 20mL, 30mL, 50mL	Power supply specification	Input voltage: 100-240 VAC
Clearance/push rate	10ml: 0.1-300ml/h , 20ml: 0.1-600ml/h	Input frequency:50/60Hz	
1/455	30ml: 0.1-900ml/h ,50ml: 0.1-1500ml/h	Battery type	11.1V 2000 mAh
Infusion mode	Speed mode,weight mode,time mode	Charging time	The charging time in the off state is more than 10h.
Speed range	0.1-1500ml/h	Display size	3.2 inches
Velocity increment	0.1ml/h	Product size	336*120*140mm
Alarm display	Audible and visual alarm	Product weight	2.0KG
Working condition	temperature:+5°C~+40°C		
	humidness:30%~80%		





# Infusion Lines - with rotary joint

"Flexible. Secure. Uninterrupted."



Prevents line twisting and ensures uninterrupted flow during transfusion—ideal for active or restless animals.



## **Subclavian Catheters**

"Precise. Secure. Specialized."

### Sizes:

Large Dog/Adult Cat: 20cm, 1.4mm Small Dog/Kitten: 18cm, 1.0mm Puppy/Neonatal: 15cm, 0.6mm

#### Features:

PE catheter & PA guidewire Sterile, single-use Radiopaque for X-ray

For central venous access



## For central venous access

### Use:

Critical care fluid therapy & medication in dogs and cats.





# Venous Butterfly Needle

"Precise. Secure. Specialized."



			The second second						
Colour	Orange 🔸	Purple	Blue •	Black	Green •	Yellow -	Cream	Pink =	
Needle Spec. (mm)	0.5×20 (256)	0.55×20 (246)	0.6×25 (230)	0.7×25 (226)	0.8×27 (216)	0.9×28 [203]	1.1×28 [196]	1.2×28 (18G)	-



C	Spec.	Oty/Pk	Oty/	ts (i)
220	Single-wing	100	5000	Non-defachable adapter
220	Single-wing	100	5000	Detachable adapter
Dout	ole-wing, Safe type	100	2000	Non-detachable adapter with safety lock protection



	Spec. Oty/Pk	Qty/Cs	(i)
20G, Pen	type 100	5000	Yellow
21G, Pen	type 100	5000	Green
22G, Pen	type 100	5000	Black
23G, Pen	type 100	5000	Blue





# **Blood Transfusion Set**

"Safe, Smooth, Reliable,"

## Length:

1.5m or 3m

## Drip Chamber:

With net filter

#### Connection:

Universal Luer lock (fits needles, catheters, IV lines)

## Material:

Medical-grade PVC & silicone (non-toxic, safe)

### Features:

Pressure-resistant, leak-proof, ready-to-use

### IV Blood Set:

Air-vented needle + 89-drop chamber + blood filter + 62PF regulator + 280cm tube + Y-port + 20cm tube + Luer lock & cap





# Blood Bags Small-Volume

"Sterile. Balanced. Ready-to-Use."



### Sizes:

50mL, 100mL, 200mL, 300mL (for small to medium animals)

## Includes:

5 bags + separate 50mL CPDA anticoagulant

## Dosing:

Pre-measured ratio (e.g., 50mL blood: 10mL anticoagulant)

#### Features:

Luer lock, medical PVC, double/single bag options

## Use:

Gently shake anticoagulant; connect via secure port





# **Vet Wound Protector**

"Secure, Sterile, Efficient,"

## Purpose:

Designed for endoscopic and small incision surgeries to retract and protect the wound.

## Key Benefits:

Optimal Surgical View - Evenly retracts the incision. Infection Prevention - Shields the wound from contamination. Procedural Efficiency - Simple and single-use.









# Blade

"Veterinary Blade-Focuses on precision for animal skin."



#### **Premium Materials**

Constructed with high-grade stainless and carbon steel for superior corrosion resistance and edge retention.

#### Multi-Sterilization Ready

Withstands repeated autoclaving, boiling, and chemical disinfection. Rust-proof treated.

## Rigorously Tested

Each blade undergoes strict quality validation including boil and performance tests. Defect-free guaranteed.

#### Trusted Globally

Widely adopted in both small and large animal surgeries for reliability and precision.





# **Animal Drainage Tube**

"Professional, Sterile, Secure,"

## **External Trauma**

6Fr 2.0×350mm 8Fr 2.6×450mm 10Fr 3.3×500mm

## **Animal Abdominal**

6Fr 2.0×350mm 8Fr 2.6×450mm 10Fr 3.3×500mm

## **Animal Thoracic**

6Fr 2.0×350mm) 8Fr 2.6×450mm) 10Fr 3.3×500mm)




# Drainage Kit - Thoracic & Abdominal

"Complete. Secure. Ready."



- 2. Guidewire
- 3. Blue-core syringe
- 4. Straight needle
- 5. Catheter fixator
- 6. Heparin cap
- 7. Dilatation catheter
- 8. Scalpel
- 9. Linker
- 10. Non-absorbable suture
- 11. Needle for injection





Procedure-Optimized – Designed for efficient thoracic/abdominal drainage.

Secure & Stable - Fixation device and suture prevent dislodgement.

Clot-Resistant - Heparin cap maintains catheter patency.

Precision Insertion - Guidewire allows accurate, minimal invasive placement.









# Silicone Drain Set

"Drain. Secure. Heal."



## I. Composition

Includes a perforated round silicone drainage tube
Includes a silicone drainage reservoir
available in three volume options: 100mL, 200mL, and 400mL
Some models come with an optional accessory: a steel trocar needle

## II. Sterilization & Packaging


Each set is individually packaged

Each package has undergone sterilization treatment

Can be directly used in aseptic operation scenarios after unpacking







## III. Material & Design

Main material: Made of biocompatible silicone

Design feature: The perforated tube allows for efficient fluid collection

## IV. Applications

Suitable for various medical scenarios

Meets the requirements of different drainage volume needs





# **Drainage Tube**

"Drain, Secure, Heal,"

#### Perforated Silicone Round

#### Size:

10FR, 12FR, 14FR, 16FR, 16FR 21 cm perforated length at the tail

#### Design:

Soft silicone, multi-hole drainage

#### Use:

Prevents fluid buildup, reduces infection

#### Features:

Flexible, biocompatible, easy to use

#### For:

Veterinary wound/surgical drainage

## Perforated Silicone Flat

#### Size:

- i. 15FR round silicone tube with 7mm × 3/4 perforated flat strip (21cm)
- ii. 15FR round silicone tube with 10mm × 3/4 perforated flat strip (21cm)

#### Cross-slit Silicone Round

#### Size:

10FR, 15FR, 19FR, 24FR Cross - slot length: 30 cm.

#### Total Length:

70-120cm

#### Use:

For efficient veterinary wound drainage

#### Features:

Soft silicone, flexible, easy to clean

#### For:

Veterinary wound/surgical drainage

#### Design:

Cross-slit opening prevents tissue ingress, allows bidirectional flow



# **Drainage Trocar**

"Complete. Secure. Efficient."







# **Jackson Pratt Drain**

"Drain. Secure. Heal."



Size(ml)	Marked Position(ml)	Remark	
100	25		1000
200	30	Clamp	
400	100		

### Material:

Medical-grade silicone (biocompatible, transparent)

## Features:

Low vacuum with anti-reflux valve Graduation markings for volume tracking Luer cap for safe and convenient use

## Compatibility:

Works with flat, fluted, or perforated drains

#### Use:

Easy-to-activate suction for efficient fluid collection





# Vet Feeding Tube (NG/OG)

"Soft, Flexible, Safe,"

Medical-grade PVC tube with radiopaque line and calibration marks for safe nasal or oral gastric intubation.



Flexible medical PVC

Sizes & Colors:

4FR (Blue), 5FR (Green), 6FR (Pink), 8FR (Clear), 10FR (Black)

## Insertion:

Nasal (NG) for conscious animals

Oral (OG) with laryngoscope for anesthetized animals (prevents tracheal misplacement)

Features: Radiopaque for X-ray verification, no stylet needed





# Vet Feeding Kit



Dual-Use Vet Feeding Kit

Designed for both pills and liquid medication.

Features an integrated syringe with soft silicone tip—secure, bite-resistant, and easy to use.





# Vet Feeding Kit - 3rd Gen

"Flexible. Kink-Resistant. Biocompatible."



2-in-1 Vet Feeding Kit

Syringe: Made of opaque white PP material

Usage: For both dry (pills) and liquid medication

Package: Two syringes per set





18cm

4.5cm

# Vet Feeding Kit - 4rd Gen

secure grip

"Flexible. Kink-Resistant. Biocompatible."

Side-Grip Design

Ergonomic side hold for easier operation

Enables simple one-handed pressing

Round handle for a comfortable and

Ergonomically Designed for Easy Use

Secure grip

One-hand press

Stable pivot

Comfortable handle

4.8cm



# **SAP Surgery Pad**

"Absorbent, Sterile, Reliable,"

## Size & Absorption:

40×60cm (>200ml)

50×60cm (>400ml)

60×90cm (≥600ml)

#### Material:

Class II medical imported cotton (23/40/50 GSM)

### Features:

Contains super absorbent polymer (SAP) 5-layer structure for optimal performance Meets most surgical needs

# **5-LAYER DESIGN**

Non-Woven Top Layer Fabric with Diamond Quilted

Super Soft Fluff Pulp Filling

Highly-Absorbent Polymer

Super Soft Fluff Pulp Filling

High-Quality Polypropylene Film Blue Backing







# Aperture Drape - Pet Sterile

"Absorbent, Sterile, Reliable,"

#### Material:

Class II medical cotton, 60 GSM hydrophilic fabric with super-absorbent polymer (SAP)


#### Design:

Adhesive-edged aperture for precise surgical access

### Sizes & Aperture Dimensions:

XS: Cover 36×43 cm | Hole 4 ×1.5 cm S: Cover 36×43 cm | Hole 9 ×3.5 cm M: Cover 36×43 cm | Hole 14 × 5 cm L: Cover 43×74 cm | Hole 20 × 7 cm XL: Cover 43×74 cm | Hole 24 × 8 cm

XXL: Cover 86×74 cm | Hole 29 ×10 cm





#### Features:

Latex-free, fluid-resistant, maintains sterile field









# Medical Gloves - PF Nitrile Exam

"Nitrile, Powder-Free, Durable,"

S (<8cm): Small hands — ideal for feline/small animal procedures

M (8-9cm): Medium hands - versatile fit for general use

L (9-10cm): Large hands — most common among veterinarians

XL (11-12cm): Extra-large hands — ensures comfort and dexterity

#### Material:

Nitrile — latex-free, puncture-resistant

#### Use:

Ideal for routine exams, surgery prep, and wound care

#### Features:

Textured grip, snug fit, safe for animal handling

#### Wear-resistant and not easy to break

Puncture-resistant/durable/not easy to be scratched during us



#### Excellent elasticity, not easy to break

Nitrile gloves have good elasticity, are strong and puncture-resistant, and fit comfortably.







# Medical Gloves - PF Latex Surgical

"Latex. Non-Sterile. Reliable."

#### Size:

S (<8cm): Small hands — ideal for feline/small animal procedures

M (8-9cm): Medium hands — versatile fit for general use

L (9-10cm): Large hands — most common among veterinarians

XL (11-12cm): Extra-large hands — ensures comfort and dexterity



#### Color:

Natural latex white only

### Packaging:

50 individually sealed pairs/box(Sterile&Non-sterile)





# Medical Gloves - PF Latex Surgical

"Pair-Packed, Sterile, Secure,"

#### Size:

S (<8cm): Small hands — ideal for feline/small animal procedures

M (8-9cm): Medium hands — versatile fit for general use

L (9-10cm): Large hands — most common among veterinarians

XL (11–12cm): Extra-large hands — ensures comfort and dexterity



#### Color:

Natural latex white only



#### Types:

Sterile (EO gas sterilized, 1 pair/bag, pre-rolled cuff)

Non-sterile (bulk packed)

#### Packaging:

50 individually sealed pairs/box(Sterile&Non-sterile)





# Positioning Mat - Vacuum

"Stable. Adaptive. Sterile."

### **Vacuum Positioning Mat Product Features**









- Exterior made of polyurethane material, durable and easy to clean.
- Radiation permeable, suitable for ultrasound, DR, CT and MRI, surgery, etc.
- Perfectly adapts to the pet's position on demand, especially for non-routine surgeries such as ophthalmology.
- High-quality safety valve that pumps in 10 seconds and stays on for over 72 hours.



#### Application Scenarios /



#### Surgery

Positioning mat can be used in surgery to fixate the surgical site of the animal, to support the surgical area, to provide a comfortable posture, to reduce surgical trauma, and to immobilize the fracture of the animal.



#### For inspection

Vacuum Positioning Mat can be used for transillumination and are suitable for ultrasound, DR, CT and MRI scenarios, making it easy to examine and analyze the animal's body condition.



### ÷

#### For rehabilitation

Vacuum Positioning Mats are often used in rehabilitation to support and stabilize the body of patients with paralysis, limb dysfunction, etc., and to help rehabilitate animals.







### Vacuum Positioning Mat Instructions for Use

- Check the positioning mat before use to ensure they are in good condition.
- Place the positioning mat horizontally on a flat surface and gently press it with your hand to make the particles evenly distributed.
- Otheck whether the pump is well connected to the connecting tube; whether the connecting tube is well connected to the air nozzle; and whether the air nozzle end of the connecting tube is connected to the air valve on the positioning mat.
- Use an air pump to pump out a small amount of air, so that the positioning mat is in a slightly tighter state, in order not to affect the positioning and shaping operations, whichever.
- The animal is laid flat on the positioning mat, the positioning mats are wrapped around the treatment area and then shaped.
- 6 After the initial molding is completed, the air pump is used to pump out part of the air in the positioning mat to form a model that basically fits the animal's body.
- According to the actual situation of the molded model to do further fine-tuning, after the adjustment is completed, use the pump again, so that the pad to achieve a vacuum state, the formation of the final model. Usually it is recommended to pump and shape in three to five times during the molding operation, which can make the product more comfortable. Reduce the difficulty of operation and better fixing effect.
- 4 After finishing shaping, remove the pump and cover the valve cover.
- After use, open the valve cover to lift the vacuum inside the pad and clean the surface.





## Precautions for use

- Please tighten the valve when pumping to avoid air leakage.
- For unanesthetized pets (cats), please work with your pet to complete nail trimming to avoid puncturing the surface of the fixation pad.
- After use, please deflate the cushion to restore the original state to reduce unnecessary tension and prolong the service life.







## **Vet Retractor**

"Retract. Expose. Secure."




#### Material:

Advanced medical polymer (non-toxic, hypoallergenic

#### Function:

Retracts wounds clearly for better visibility and shorter surgery time

## **Hook Types:**







# **SUB Implantable Consumable**

"Innovative. Indwelling. Integral."

#### Intended Use

The Vetmediux Drug Support Port provides reliable subcutaneous access for long-term drug therapy and fluid management in companion animals. It enables repeated infusion and sampling while eliminating the need for repeated venipuncture.

#### Core Components

- 1. Implantable Port: Features a self-sealing silicone septum within a biocompatible housing
- 2. Specialized Catheter: 1.8mm × 750mm catheter with radiopaque markers
- 3. Procedure Kit: Complete set for implantation and maintenance procedures

#### Key Advantages

- 1. Long-Term Solution: Designed for extended indwelling periods with minimal complications
- 2. Enhanced Safety: Closed system reduces infection risks and simplifies access
- 3. Universal Compatibility: Available configurations suit patients of varying sizes

#### Maintenance Protocol

- 1. Follow flushing schedule: pre-discharge, 1 week, 1 month, then quarterly
- 2. Utilize dedicated flushing solutions and techniques
- 3. No scheduled replacement when functioning properly

#### Clinical Monitoring

The Vetmediux Drug Support Port remains in place indefinitely while maintaining proper function. Regular assessment during quarterly flushing sessions determines continuing suitability. Replacement decisions integrate individual patient factors and operational status.

#### Artical number:



1) VLSUB01



2) VLSUB02



3) VLSUB03



4) VLSUB04



5) VLSUB05

- Veterinary Implantable Drug Delivery Port 1.8×750mm
- Veterinary Subcutaneous Ureteral Bypass (SUB) Unilateral Type A
- Veterinary Subcutaneous Ureteral Bypass (SUB) Unilateral Type B
- Veterinary Subcutaneous Ureteral Bypass (SUB) Bilateral Type C
- Veterinary Ureteral Flushing Kit

A veterinary-specific implantable device.





# Drug Delivery Port - 1.8×750mm

"Innovative. Indwelling. Integral."

### Accessories:

- 1. 1x TC Titanium Alloy
   Injection Port Body
- 2. 1x Intravenous Infusion Tube
- 1x 2.0×231mm Introducer

  Needle
- 4. 1x Blue Protective Cover
- 5. 1x Alloy Guide Wire
- 1x Splittable Catheter Sheath
- 7. 1x Heparin Cap
- 2x 0.7\*32mm TWLB
   Injection Needles, 1x 18G
   Injection Needle
- 9. 2x 10ml Disposable Sterile Syringes
- 10. 1x 18G×6.5cm Straight
   Puncture Needle
- 11. 1x Disposable
   Plastic-Handle Scalpel
- 12. 1x Vessel Dilator
- 13. 1x 20G Disposable
   Implantable Port Access Needle








# SUB Ureteral Bypass - Unilateral Type A

"Innovative. Indwelling. Integral."

### Accessories:

- 1. 1x 7F × 231mm Bladder
   Straight Catheter + Hollow Rigid
   Tube
- 2. 1x 6.5F × 200mm J-Tube with Fixation Ring + Hollow Rigid Tube
- 1x TC Titanium Alloy Straight Connector
- 4. 2x Blue Protective Covers
- 5. 1x Alloy Guide Wire
- 6. 2x 22G × 0.75 Huber Needles
- 7. 1x Heparin Cap
- 2x 18G Disposable Floating Needles
- 9. 1x 0.25ml Tissue Adhesive
- 10. 1x Glue Dropper
- 11. 1x Disposable
   Plastic-Handle Scalpel
- 12. 1x Vessel Dilator







# SUB Ureteral Bypass - Unilateral Type B

"Innovative. Indwelling. Integral."

#### Accessories:

- 1x TC Titanium Alloy Injection Port Body
- 1x 6.5F × 200mm J-Tube with Fixation Ring + Hollow Rigid Tube
- 1x 7F × 231mm Bladder
   Straight Catheter + Hollow Rigid
   Tube
- 1x TC Titanium Alloy 3-Way
  Connector
- 5. 5x Blue Protective Covers
- 6. 1x Alloy Guide Wire
- 7. 1x Port Connection Tube
- 8. 2x 22G × 0.75 Huber

#### Needles

- 9. 1x Heparin Cap
- 2x 18G Disposable Floating Needles
- 11. 1x 0.25ml Tissue Adhesive
- 12. 1x Adhesive Applicator Dropper
- 13. 1x Disposable
- Plastic-Handle Scalpel
- 14. 1x Vessel Dilator
- 15. 1x 20G Disposable Implantable Port Access Needle







# SUB Ureteral Bypass - Unilateral Type C

"Innovative. Indwelling. Integral."

#### Accessories:

- 1x TC Titanium Alloy Implantable Port Body
- 1x 7F × 231mm J-Tube with Fixation Ring + Hollow Rigid Tube
- 1x 7F × 231mm J-Tube with Fixation Ring + Hollow Rigid Tube
- 1x 6.5F × 200mm Bladder Straight Catheter + Hollow Rigid Tube
- 1x TC Titanium Alloy 4-Way Connector
- 6. 6x Blue Protective Covers
- 7. 1x Alloy Guide Wire
- 8. 1x Port Connection Tube
- 2x 22G × 0.75 Huber Needles
- 10. 1x Heparin Cap
- 2x 18G Disposable Floating Needles
- 12. 1x 0.25ml Tissue Adhesive
- 13. 1x Adhesive Applicator Dropper
- 14. 1x Disposable
   Plastic-Handle Scalpel
- 15. 1x Vessel Dilator
- 16. 1x 20G Disposable Implantable Port Access Needle








# **Ureteral Flushing Kit**

"Innovative. Indwelling. Integral."



\*(Package also includes 2% EDTA irrigation solution)\*

## Accessories:

- 1. 1x Three-Way Extension Tube
- 2. 1x 10ml Disposable Sterile Luer Lock Syringe
- 3. 1x 5ml Disposable Sterile Syringe
- 4. 1x 20G Disposable Implantable Port Access Needle



# Injection Needles - Disposable

"Sharp. Sterile. Secure."

ŘmmxLmm	Gauge	Colour
0.5 x 16 mm	25G x 5/8"	Orange
0.6 x 16 mm	25G x 5/8"	Light Blue
0.6 x 25 mm	,23G x 1"	Blue
0.8 x 25 mm	21G x 1"	Light Green
0.8 x 10 mm	21G x 2/1"	Green
0.8 x 40 mm	21G x 12/1°	Teal
0.8 x 16 mm	25G x 5/8"	Yellow
0.9 x 25 mm	20G x 1"	Light Yellow
0.9 x 40 mm	20G x 12/1"	Beige
1.1 x 25 mm	19G x 1"	Transparent
1.1 x 40 mm	19G x 12/1"	Transparent
1.1 x 16 mm	19G x 5/8"	Transparent
1.2 x 40 mm	19G x 11/4"	Pink"
1.25 x 25 mm	18G x 1"	Light Pink
1.6 x 25 mm	16G x 1"	Light Peach
1.6 x 40 mm	16G x 12/1*	White (empty circle)
0.6 x 30 mm	23G x 12/1"	Gray
0.7 x 30 mm	22G x 19/1	Dark Gray
0.4 x 20 mm	27G x 4/5"	Brown
0.45 x 25 mm	28G x 1"	Dark Brown
0.55 x 25 mm	24G x 1"	Purple
0.7 x 20 mm	22G x 2/1"	Lime
0.9 x 13 mm	20G x 2/1"	Gold
0.5 x 10 mm	25G x 3/8"	Dark Brown
0.6 x 19 mm	23G x 3/8"	Light Blue

## Design:

High-quality stainless steel, ultra-sharp facet point

#### Benefits:

Smooth puncture, minimal tissue damage, reduced pain

#### Hub:

Semi-transparent for easy reflux observation

## Sterilization:

EO gas sterilized, paper-film peel packaging

## Size Range:

Comprehensive options for veterinary use

## Packaging:

100 needles per pack







## Hypodermic Needles - 304 Stainless Steel

"Sharp. Sterile. Secure."



Applicable Objects	Select Needle
	NO.7
<b>₩</b>	NO.7 \ NO.9
77	NO.7 \ NO.9 \ NO.12 \ NO.16
APP	NO.12 \ NO.14 \ NO.16
Experiment	NO.7 \ NO.9 \ NO.12(More commonly used)
Mechanical filling with oiling agent	NO.14 \ NO.16

Material: Premium stainless steel (14G-22G)

Use: Reliable for vaccines and medications in livestock

#### Features:

Elongated design for deep penetration Thickened tube resists bending Reinforced needle seat prevents loosening



# Sterile Syringe - Needle-Free

"Accurate. Smooth. Versatile."







#### Luer Lock

#### Luer Slip

The fip's sleek surface ensures comfortable direct use.
Attachments are made effortlessly fhrough a simple sliding mechanism.
Optional capacities: I.M., Z.W., J.M., S.W., I.D.M., Z.O.M.

## Large Slip Tip

For measuring and dispensing fluids/chemicals, lab use & experiments, succulents garden watering. Optional capacities: 20ML, 60ML, 100ML, 150ML, 200ML, 300ML,350ML,60ML,

Specification	Whole length (mm)	Tube length (mm)	Tube diameter (mm)
1ml	80	55	12
2.5ml	85	65	13
3ml	85	65	13
5ml	85	58	15
10ml	100	75	20
20ml	120	90	25
30ml	135	105	25
60ml	145	110	35



### Sizes:

1ml, 2ml, 3ml, 5ml, 10ml, 20ml, 50ml

#### Features:

Single-use, sterile-packed

Precise dosing scale

Smooth plunger, easy to use

Durable and lightweight












## Fixed-Needle Syringe - Luer Slip / Luer Lock

"Ready. Precise. Efficient."



MATERIAL	Plastic and Metal
CAPACITY	1ml, 3ml, 5ml, 10ml, 20ml, 30ml, 50ml
PACKAGE	Box or Bag
FEATURE	Disposable
USEFOR	Veterinary, animal, poultry, food





## Luer slip adapter:

Description	Size	Graduation	ømm x Lmm	Gauge	Colour
Insulin Syringe with fixed needle, centric (USO)	0.5ml		0.33 x 12 mm	29Gx/4"	White
Insulin Syringe with fixed needle, centric (U100)	1ml		0.33 x 12 mm	29Gx%"	White
TBC Syringe with fixed needle, centric	1ml		0.5 x 16 mm	25G x %"	Orange
Syringe with fixed needle, centric	1ml		0.5 x 16 mm	25G x %"	Orange
Syringe with fixed needle, centric	1ml		0.45 x 40 mm	26G x 115"	Brown
Syringe with fixed needle, centric	2 ml		0.8 X 12 mm	21G x 1/2"	Green
Syringe with fixed needle, centric	2 ml		0.6 x 25 mm	33 x 1"	Light Blue
Syringe with fixed needle, centric	2 ml		0.7 x 32 mm	22G x 1%"	Black
Syringe with fixed needle, centric	2 ml		0.8 x 40 mm	210 x 1½"	Dark Green
Syringe with fixed needle, centric	2 ml	3 mi	0.8 x 16 mm	21G x %"	Light Green
Syringe with fixed needle, centric	5 ml	6 ml	0.7 x 32 mm	22G x 114"	Black
Syringe with fixed needle, centric	5 ml	6 ml	0.8 x 40 mm	21G x 1½"	Light Green
Syringe with fixed needle, centric	5 ml	6 mi	125 x 40 mm	18G x 1%"	Pink
Syringe with fixed needle, eccentric	10 ml	12 ml	0.8 x 40 mm	21G x 1½"	Light Green
Syringe with fixed needle, eccentric	10 ml	12 ml	0.9 x 40 mm	20G x 11/4"	Yellow
Syringe with fixed needle, eccentric	10 ml	12 ml	0.7 x 40 mm	22G x 115"	Black
Syringe with fixed needle, eccentric	20 ml	25 ml	1.25 x 40 m	18G x 1½"	Light Pink

## Luer lock adapter:

Description	Size	Graduation	ømm x Lmm	Gauge	Colour
Syringe with fixed needle, centric	1 ml		0,5 x 16 mm	25G x %"	Orange
Syringe with fixed needle, centric	2 ml	3 ml	0.7 x 32 mm	22G x T¼"	Black
Syringe with fixed needle, centric	2 ml	3 ml	.8 x 40 mm	21G x 114"	Green
Syringe with fixed needle, centric	5 mll	6 mll	2 x 40 mm	18G x 1½"	Pink








# Insulin Syringes - Disposable

"safety. accuracy. animal comfort."



Ultra-Sharp Needles: Silicon-coated for smooth insertion and reduced pain. Clear Barrel Design: High transparency with precise scale markings for easy and accurate dosing.

## Complete Range:

From standard syringes (1mL–50mL) to specialty insulin and TB syringes, all available in various needle gauges (e.g. 16G–26G).

Our veterinary syringes and needles are engineered for safety, accuracy, and animal comfort. Made from medical-grade PP material, each product is disposable, sterile, and fully complies with CE and ISO standards.





# 1mL TB Syringe - Optional Needle

"safety. accuracy. animal comfort."



### **Dual Connector:**

Available with both Luer Lock (secure twist-on) and Luer Slip (smooth push-on) options.

## Precise Dosing:

Clear scale markings enable accurate medication measurement.

### Leak-Resistant:

Smooth plunger prevents fluid leakage.

## Sterile & Single-Use:

Individually wrapped and capped for safety.

Purpose: For veterinary use only.







# Blood Gas Syringe - Arterial

"Precise, Critical, Reliable,"



### Purpose:

Critical for diagnosing respiratory, circulatory, and metabolic diseases by measuring blood oxygen, acid-base balance, and electrolyte levels. Essential for surgery monitoring and emergency care.

## Key Features:

Precise Anticoagulant: Calcium-balanced lithium heparin (30 IU/mL for 1mL; 80 IU/mL for 3mL) prevents clotting without distorting calcium test results.

## Common Specs:


1mL & 3mL sizes
With/without needle (23G×25mm if equipped)

#### Clinical Value:

Provides the only reliable data for detecting hypoxemia and acid-base imbalance.







Use calcium-balanced heparin lithium as anticoagulant Blood collection with default value, very convenient to operate New medical device in arterial blood specimen collection

#### Sharper needle point

The special grinding and labrication technology can ease the pain from the patients.

#### Self-venting Function

This function can make sure that the air inside the syringer being exhaused completely after arteral blood entering thems strend markets.

#### Hermetic type

The internal blood will block the stotem when petting in touch with the justion. Wevent the air from petting

#### Blood collection needle sealing

Optimized butyl rubber can seal fre needle and prevent needle stick injuries.



#### Tranparent needle sea

Transparent needle seat can help the manes to observe the blood return clearly and make accurate judgment on puncture.

#### Sprayed, solid calciumbalanced heparin lithiun

The solid calcum-halanced heparin lithium is sprayed equiliby inside the syringe, which can make sure of the accuracy of the detected.

#### Conchead sealing cover

The double-by et design can seal the somple completely and polyte the ar-


#### Simplify the operation

Save W7% operation time and simplify the





# Vet Irrigation Syringe



## Precision Tips:

- 0.45mm/16mm/50°
- 0.5mm/32mm/30°
- 0.7mm/32mm/30°

## Broad-Spectrum Use:

Ideal for delicate flushing of eyes, tear ducts, teeth, ears, and nose

## Ready-to-Use:

EO sterilized and individually packaged

## Ergonomic Design:

Angled tip enables comfortable and accurate application





# **Urine Test Strips**


"Accurate, Certified, Stable"


Accuracy: 96%
Shelf Life: 2 years
Certifications: CE, ISO
Material: PET

Language : English (customization available)

Test Time : 60 seconds

Advantage : Resistant to VC (ascorbic acid) interference





### **Test Parameters:**

Leukocytes, Nitrite, Urobilinogen, Protein, pH, Occult Blood, Specific Gravity, Ascorbic Acid, Ketones, Bilirubin, Glucose. Intended Use: For urine analysis.







# Test cards

Parameter A	bbreviation	Principle	Reference Range
PH	PH	Acid-base indicator method	PH4.5-8.0
Specific Gravity	SG	Polyelectrolyte ion dissociation method	1.015-1.025
Protein	PRO	PH indicator protein error method	Negative
Glucose	GLU	Glucose oxidase peroxidase method	Negative
Bilirubin	BIL	Azo reaction method	Negative
Urinary bile proto	URO	Aldehyde reaction, diazotization method	Negative
Ketone	KET	Sodium nitroso ferricyanide method	Negative
Nitrite	NIT	Nitrite reduction	Negative
Blood or Red Blood Cells	BLD	Hemoglobin peroxidase method	Negative
White Blood Cell	LEU	Esterase method	Negative
Vitamin C	VitC	Indole enzyme method	Negative
*Urinary Creatinine	Cr	3,5-two Nitrobenzoic Acid Method	Negative
*Urinary calcium	Ca	Phthalo phthalein complex ketone method	Negative
*Urinary Microalbuminuria	uMA	Fluorescein dye method	Negative





#### SUMMARY

Urinalysis Reagent Strips are made for urinalysis of both qualitative and semi-quantitative, which are in vitro reagent for diagnostics. It tests Leukocytes, Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ascorbic Acid, Ketone, Bilirubin, Glucose, Microalbumin, Creatinine, Calcium in urine. Please refer to the out-side box and bottle label for the specific test parameters of the product you are using. Please read this direction carefully before using.

The results on the strips can be read visually and instrumentally

### **SPECIFICATIONS**

Model No.	Box size	Pcs/box	N.W. (kg)	G.W. (kg)
URS-14	53*53*147mm	100	0.12	0.15

#### SPECIMEN COLLECTION AND PREPARATION

Use only clean dry container to collect urine and should be shocked before testing and test it within 2 hours. Any operations must be in the clean environment.

#### **TEST CONDITIONS**

Ambient temperature: 20C-30C, relative humidity<80%, the best test temperature: 23 C -27 C .

#### **STORAGE**

Store between 2-30 tin dry condition. Keep away from refrigerator direct sunlight Do not touch test area of reagent strips. Isolated from damp, light and high temperature for the aim of preserving the reaction activity of reagent.





#### TEST PROCEDURE

- 1. Remove one strip from the bottle and replace the cap immediately.
- Immerse the reagent area of the strip in the urine specimen and take it out quickly.
- 3. Wipe off excess urine against the rim of the specimen container.
- 4.Read the test results carefully within 60 seconds in a good light and with the test area held near the appropriate color chart on the bottle label. Changes in color that appear only along the edges of the test pads or after moving than 2 minutes have passed are of no diagnostic significance. Results with leukocytes test portion can be read within 120 seconds.

If reading instrumentally, carefully follow the directions given in the appropriate instrument operating manual.

#### ATTENTION

Water cannot be used as negative quality control liquid ,Antiseptic of urine cannot prevent the ketone, bilirubin and urobilinogen from deteriorated. For the long time urine specimen, the test results of glucose, pH, nitrite and blood can be affected cause of bacterial growth.







#### WARNINGS AND PRECAUTIONS

- 1.Do not remove desiccant from the bottle.
- Do not touch test area of Urine Reagent Strips.
- 3.Do not open container until ready to use.
- 4.The use of urine preservatives can prevent the decomposition of ketone, bilirubin and urobilinogen in the urine.
- 5.Do not store the sample long time (one hour or longer) before testing.

#### LIMITATION

Comparison to the color chart is dependent on the interpretation of the individual. It is therefore, recommended that all laboratory personnel interpreting the results of these strips be tested for color blindness. As with all laboratory tests, definitive diagnostic or therapeutic decisions should not be based on any single test or method.



# **Animal Urea Analyzer**

"Sterile, Versatile, Precision"



## FREEZING REFRACTOMETER


Vehicle urea concentration: 30 ~ 35%

Specific gravity of battery liquid: 1.1 ~ 1.4

Freezing point of antifreeze: - 50 ~ 0 °C

Glass cleaning agent: - 40 ~ 0 °C











#### Plating and soldering process


High hardness and wear resistance, anti-corrosion and anti rust, bright as new, atmospheric and beautiful



Soft feel, beautiful and comfortable, light rotation and not easy to slide

#### Clear reading

Fast detection, accurate measurement, clear display and convenient use



## High definition prism



## Calibration screw











# **Biopsy Punch**

"Sterile. Sharp. Precise."



00G(10mm) 0G(8mm) 2G(6mm) 4G(5mm) 6G(4mm) 8G(3mm) 12G(2mm) 14G(1.5mm)

## Precision Sizes:

1.5-10mm (1.5/2/3/4/5/6/7/8/9/10mm)

## EO Sterilized:

Individually packed, single-use

## Sharp & Reliable:

Stainless steel, suitable for skin procedures and lab use

The operation is simple









# Rapid Tests

"Fast. Accurate. Essential."

# **!** Fenine Only

No.	Abbr.	Description	Sample
НАЗ8	FPV Ag	Feline Panleucopenia Antigen Rapid Test	Feces
HA39	FELV Ag	Feline Leukemia Virus Antigen Rapid Test	Whole/Serum/Plasma
HA40	FIV Ab	Feline Immuno deficiency Virus Antibody Rapid Test	Whole/Serum/Plasma
HA41	FIV Ab-FELV Ag	FIV Ab/FeLV Ab Combo Rapid Test	Whole/Serum/Plasma
HA42	FCOV Ab	Feline Coronavirus Antibody Rapid Test	Whole/Serum/Plasma
HA43	FCV Ag	Feline Calicivirus Antigen Rapid Test	Oral, nasal, and eye secretions
*HA44	FRL	Feline Pregnancy Relaxin Test (FRL)	Whole/Serum/Plasma
HA45	FIPV Ab	Feline Infectious Peritonitis Virus Antibody Rapid Test	Serum/Plasma/Ascites
HA46	FCoV Ag	Feline Coronavirus Antigen Rapid Test	Feces
HA47	FIV-FELV-HW	FIV Ab/Felv Ag/Heartworm Antigen Combo Rapid Test	Whole/Serum/Plasma
HA48	FPV-FCoV-GIA	Feline Panleukopenia/ Corona/Giardia Antigen Combo Rapid Test	Feces
HA49	FHV Ag	Feline Herpesvirus Antigen Rapid Test	Oral, nasal, and eye secretions
HA50	FCV-FHV(Type-1)Ag	Feline Calicivirus/Herpesvirus Type-1 Antigen Rapid Test	Oral, nasal, and eye secretions
HA51	FCV Ab	Feline Caficivirus Antibody Rapid Test	Whole/Serum/Plasma
HA52	FPV Ab	Feline Panleukopenia Virus Antibody Rapid Test	Whole/Serum/Plasma
HA53	FHV Ab(Type-1)	Feline Herpesvirus Type-1 Antibody Rapid Test	Whole/Serum/Plasma
HA54	FIV-FELV-FCOV-FHW	FIV Ab/Felv Ag/FCoV Ag/Heartworm Ag Combo Rapid Test	Whole/Serum/Plasma
HA55	FPV-FCV-FCOV-GIA	Feline Panleukopenia/ Calicivirus/Corona/Giardia Combo Rapid Test	Oral, nasal, eye secretions and feces
HA56	FHV-FCV-FPV Ab	FHV/FCV/FPV Antibody Combo Rapid Test	Whole/Serum/Plasma
HA57	FELV Ag(Sallvia)	Feline Leukemia Virus Antigen Rapid Test(Salvia)	Salvia
HA58	FIV Ab(Salvia)	Feline Immuno deficiency Virus Antibody Rapid Test(Salvia)	Salvia
HA59	FIV Ab-FELV Ag(Salvia)	FIV Ab/FeLV Ag Combo Rapid Test(Salvia)	Salvia
HA60	FMP	Feline Mycoplasma Rapid Test	Oral, nasal, and eye secretions
HA61	FPV-FCoV Ag	Feline Panleukopenia/ Corona Combo Rapid Test	Feces
HA62	F.Blood typing	Feline Blood typing Card	Whole Blood
HA76	FIPV Ag	Feline Infectious Peritonitis Virus Antigen Rapid Test	Whole Blood/Serum/Plasma
HA85	FIPV Ag	Feline Infectious Peritonitis Virus Antigen Rapid Test	feces, rectal swab/ascites
HA86	FHV-FCV-FMP Ag	FHV-FCV-FMP Ag Combo Rapid Test	Oral, nasal, and eye secretions

# Pet Rpid Test



Abbr.	Description	Sample
CDV Ag	Canine Distemper Virus Antigen Rapid Test	Oral, nasal, and eye secretions
CPV Ag	Canine Parvo Virus Antigen Rapid Test	Feces
CHW Ag	Canine Heartworm Antigen Rapid Test	Whole/Serum/Plasma
CCV Ag	Canine Coronavirus Antigen Rapid Test	Feces
CPV Ab	Canine Parvovirus Antibody Rapid Test	Whole/Serum/Plasma
CDV Ab	Canine Distemper Virus Antibody Rapid Test	Whole/Serum/Plasma
C.Bru Ab	Canine Brucella Antibody Rapid Test	Whole/Serum/Plasma
CPV-CCV Ag	Canine Parvo/Corona Antigen Combo Rapid Test	Feces
ANA Ab	Canine Anaplamsa Antibody Rapid Test	Whole/Serum/Plasma
CPV-CCV-Giardia Ag	CPV/CCV/Giardia Antigen Combo Rapid Test	Feces
LSH Ab	Canine Leishmania Antibody Rapid Test	Whole/Serum/Plasma
EHR-LYM-ANA Ab	Ehrlichia /LymeiAnoplasma Antibody Combo Rapid Test	Whole/Serum/Plasma
Lyme Ab	Canine Lyme Disease Antibody Rapid Test	Whole/Serum/Plasma
CRV Ag	Canine Rotavarus Antigen Rapid Test	Feces
CPV-CCV-CRV Ag	CPV/CCV/CRV Antigen Combo Rapid Test	Feces
CDV-CPV Ag	Canine Distemper/Parvo Virus Antigen Combo Rapid Test	Feces
CDV-CAV Ag	Canine Distemper/Adenovirus Antigen Combo Rapid Test	Oral, nasal, and eye secretion
CIV Ag	Canine Influenza Virus Antigen Rapid Test	Oral, nasal, and eye secretion
CDV-CAV-CIV Ag	Canine Distemper/Adeno/Influenza Combo Rapid Test	Oral, nasal, and eye secretion
RLN	Canine Pregnancy Relaxin Rapid Test	Whole/Serum/Plasma
CPV-CDV Ab	CPV/CDV Antibody Combo Rapid Test	Whole/Serum/Plasma
EHR-ANA Ab	Ehrlichia/Anaplasma Antibody Combo Rapid Test	Whole/Serum/Plasma
BG Ab	Babesia Gibsoni Antibody Rapid Test	Whole/Serum/Plasma
EHR Ab	The second secon	
	Ehrichia Antibody Rapid Test	Whole/Serum/Plasma
EHR-BG-ANA Ab	Ehrlichia/Babesia Gibsoni/Anaplasma Antibody Combo Rapid Test	Whole/Serum/Plasma
EHR-LYM-ANA-CHW	Ehrlichia/Lyme/Anaplasma/Heartworm Combo Rapid Test	Whole/Serum/Plasma
EHR-BG-ANA-CHW Ab	Ehrlichia/Babesia Gibsoni/Anaplasma/Heartworm Combo Rapid Test	Whole/Serum/Plasma
CPV-CCV-GIA-CRV Ag	CPV/CCV/GIA/CRV Combo Rapid Test	Feces
BC Ab	Babesia Canis Antibody Rapid Test	Whole/Serum/Plasma
CAV Ag	Canine Adenovirus Antigen Rapid Test	Oral, nasal, and eye secretion
CPfV Ag	Canine Parainfluenza Virus Ag Rapid Test	Oral, nasal, and eye secretion
EHR-LYM-ANA-BAB Ab	Ehrlichia-Lyme-Anaplasma-Babesia Gibsoni Antibody Combo Rapid Test	Whole/Serum/Plasma
CDV-CAV-CIV-CPIV Ag	CPV/CAV/CIV/CPIV Antigen Combo Rapid Test	Oral, nasal, and eye secretion
CDV-CPV-CCV-GIA Ag	CDV/CPV/CCV/GIA Antigen Combo Rapid Test	Oral, nasal, and eye secretion
EHR-LSH Ab	Ehrlichia/Leishimania Antibody Combo Rapid Test	Whole/Serum/Plasma
EHR-ANA-BG-BC Ab	EHR/ANA/BG/BC Antibody Combo Rapid Test	Whole/Serum/Plasma
cProg	Canine Progestrone Rapid Test	Serum/Plasma
EHR+LSH+ANA+CHW	EHR Ab+LSH Ab+ANA Ab+CHW Ag Combo Rapid Test	Whole Blood/Serum/Plasma
C.Bru Ag	Canine Brucella Antigen Rapid Test	Semen, vaginal secretions
CAV+CDV+CPV Ab	CAV+CDV+CPV Ab Combo Rapid Test	Whole Blood/Serum/Plasma
Angiostrongylus vasorum Ag	Angiostrongylus vasorum Ag Rapid Test	Whole Blood/Serum/Plasma
FOB	Fecal Occult Blood Test	Feces
Anaplasma Ab+Dirofilana Ag	Anaplasma Ab+Diroffaria Ag+Ehrlichia Ab+Borrella Ab Combo Rapid Test	Whole Blood/C Dis-
+Ehrlichia Ab+Borrelia Ab	remplaced Performance agreement Auroportona Auroporto Rapid Test	viriole biood/serum/Plasma
Cytomegalovirus Ag	Commence to Commence to	
+Parvovirus Ag	Cytomegatevirus Ag+Parvevirus Ag	
CPV-CCV-GIA-CRY Ag	Giardia+PARVO Vitrus-Cryposporidium+Corona virus Combo Rapid Test	Feces
EHR-ANA-Lyme-BAB-LSH	Ehlrichia+Anaplasma+Lyme disease+Babesia+Leishmania Combo Rapid Test	Whole Blood/Serum/Plasma

# Pet Rpid Test




Abbr.	Description	Sample
CDV Ag	Canine Distemper Virus Antigen Rapid Test	Oral, nasal, and eye secretions
CPV Ag	Canine Parvo Virus Antigen Rapid Test	Feces
CHW Ag	Canine Heartworm Antigen Rapid Test	Whole/Serum/Plasma
CCV Ag	Canine Coronavirus Antigen Rapid Test	Feces
CPV Ab	Canine Parvovirus Antibody Rapid Test	Whole/Serum/Plasma
CDV Ab	Canine Distemper Virus Antibody Rapid Test	Whole/Serum/Plasma
C.Bru Ab	Canine Brucella Antibody Rapid Test	Whole/Serum/Plasma
CPV-CCV Ag	Canine Parvo/Corona Antigen Combo Rapid Test	Feces
ANA Ab	Canine Anaplamsa Antibody Rapid Test	Whole/Serum/Plasma
CPV-CCV-Giardia Ag	CPV/CCV/Giardia Antigen Combo Rapid Test	Feces
LSH Ab	Canine Leishmania Antibody Rapid Test	Whole/Serum/Plasma
EHR-LYM-ANA Ab	Ehrlichia /Lymn/Anaplasma Antibody Combo Rapid Test	Whole/Serum/Plasma
Lyme Ab	Canine Lyme Disease Antibody Rapid Test	Whole/Serum/Plasma
CRV Ag	Canine Rotavirus Antigen Rapid Test	Feces
CPV-CCV-CRV Ag	CPV/CCV/CRV Antigen Combo Rapid Test	Feces
CDV-CPV Ag	Canine Distemper/Paryo Virus Antigen Combo Rapid Test	Feces
CDV-CAV Ag	Canine Distemper/Adenovirus Antigen Combo Rapid Test	Oral, nasal, and eye secretions
CIV Ag	Canine Influenza Virus Antigen Rapid Test	Oral, nasal, and eye secretions
CDV-CAV-CIV Ag	Canine Distemper/Adeno/Influenza Combo Rapid Test	Oral, nasal, and eye secretions
RLN	Canine Pregnancy Relaxin Rapid Test	Whole/Serum/Plasma
CPV-CDV Ab	CPV/CDV Antibody Combo Rapid Test	Whole/Serum/Plasma
EHR-ANA Ab	Ehrichia/Anaplasma Antibody Combo Rapid Test	Whole/Serum/Plasma
BG Ab	Babesia Gibsoni Antibody Rapid Test	Whole/Serum/Plasma
EHR Ab	Ehriichia Antibody Rapid Test	Whole/Serum/Plasma
EHR-BG-ANA Ab	Ehrlichia/Babesia Gibsoni/Anaplasma Antibody Combo Rapid Test	Whole/Serum/Plasma
EHR-LYM-ANA-CHW	Ehrlichia/Lyme/Anaplasma/Heartworm Combo Rapid Test	Whole/Serum/Plasma
EHR-BG-ANA-CHW Ab	Ehrlichia/Babesia Gibsoni/Ariaplasma/Heartworm Combo Rapid Test	Whole/Serum/Plasma
CPV-CCV-GIA-CRV Ag	CPV/CCV/GIA/CRV Combo Rapid Test	Feces
BC Ab	Babesia Canis Antibody Rapid Test	Whole/Serum/Plasma
CAV Aq	Canine Adenovirus Antigen Rapid Test	Oral, nasal, and eye secretions
CPIV Ag	Canine Parainflueriza Virus Ag Rapid Test	Oral, nasal, and eye secretions
EHR-LYM-ANA-BAB Ab	Ehrlichia-Lyme-Anaplasma-Babesia Gibsoni Antibody Combo Rapid Test	Whole/Serum/Plasma
CDV-CAV-CIV-CPIV Ag	CPV/CAV/CIV/CPIV Antigen Combo Rapid Test	Oral, nasal, and eye secretion:
CDV-CPV-CCV-GIA Ag	CDV/CPV/CCV/GIA Antigen Combo Rapid Test	Oral, nasal, and eye secretion
EHR-LSH Ab	Ehrlichia/Leishimania Antibody Combo Rapid Test	Whole/Serum/Plasma
EHR-ANA-BG-BC Ab	EHR/ANA/BG/BC Antibody Combo Rapid Test	Whole/Serum/Plasma
Prog	Canine Progestrone Rapid Test	Serum/Plasma
EHR+LSH+ANA+CHW	EHR Ab+LSH Ab+ANA Ab+CHW Ag Combo Rapid Test	Whole Blood/Serum/Plasma
C,Bru Ag	Canine Brucella Antigen Rapid Test	Semen, vaginal secretions
CAV+CDV+CPV Ab	CAV+CDV+CPV Ab Combo Rapid Test	Whole Blood/Serum/Plasma
Angiastrongylus vasorum Ag	Angiostrongylus vasorum Ag Rapid Test	Whole Blood/Serum/Plasma
FOB	Fecal Occult Blood Test	Feces
Anaplasma Ab+Dirofilaria Ag		
+Ehrlichia Ab+Borrelia Ab	Anaplasma Ab+Diroftana Ag+Etrilichia Ab+Borrelia Ab Combo Rapid Test	Whole Blood/Serum/Plasma
Cylomegalovirus Ag +Parvovirus Ag	Cytomogaloveus Ag+Parvovirus Ag	
CPV-CCV-GIA-CRY Ag	Giardis+PARVO Vitrus-Cryposporidium+Corona virus Combo Rapid Test	Feces
EHR-ANA-Lyme-BAB-LSH	Ehlrichia+Anapiasma+Lyme disease+Babesia+Leishmania Combo Rapid Test	

# Pet Rpid Test



Abbr.	Description	Sample
CDV Ag	Canine Distemper Virus Antigen Rapid Test	Oral, nasal, and eye secretion
CPV Ag	Canine Parvo Virus Antigen Rapid Test	Feces
CHW Ag	Canine Heartworm Antigen Rapid Test	Whole/Serum/Plasma
CCV Ag	Canine Coronavirus Antigen Rapid Test	Feces
CPV Ab	Canine Parvovirus Antibody Rapid Test	Whole/Serum/Plasma
CDV Ab	Canine Distemper Virus Antibody Rapid Test	Whole/Serum/Plasma
C <sub>s</sub> Bru Ab	Canine Brucella Antibody Rapid Test	Whole/Serum/Plasma
CPV-CCV Ag	Canine Parvo/Corona Antigen Combo Rapid Test	Feces
ANA Ab	Canine Anaplamsa Antibody Rapid Test	Whole/Serum/Plasma
CPV-CCV-Giardia Ag	CPV/CCV/Giardia Antigen Combo Rapid Test	Feces
LSH Ab	Canine Leishmania Antibody Rapid Test	Whole/Serum/Plasma
EHR-LYM-ANA Ab	Ehrlichia /Lyme/Anaplasma Antibody Combo Rapid Test	Whole/Serum/Plasma
.yme Ab	Canine Lyme Disease Antibody Rapid Test	Whole/Serum/Plasma
CRV Ag	Canine Rotavirus Antigen Rapid Test	Feces
CPV-CCV-CRV Ag	CPV/CCV/CRV Antigen Combo Rapid Test	Feces
CDV-CPV Ag	Canine Distemper/Parvo Virus Antigen Combo Rapid Test	Feces
CDV-CAV Ag	Canine Distemper/Adenovirus Antigen Combo Rapid Test	Oral, nasal, and eye secretion
CIV Ag	Canine Influenza Virus Antigen Rapid Test	Oral, nasal, and eye secretion
CDV-CAV-CIV Ag	Canine Distemper/Adeng/Influenza Combo Rapid Test	Oral, nasal, and eye secretion
RLN	Canine Pregnancy Relaxin Rapid Test	Whole/Serum/Plasma
CPV-CDV Ab	CPV/CDV Antibody Combo Rapid Test	Whole/Serum/Plasma
EHR-ANA Ab	Ehrlichia/Anaplasma Antibody Combo Rapid Test	Whole/Serum/Plasma
3G Ab	Babesia Gibsoni Antibody Rapid Test	Whole/Serum/Plasma
EHR Ab	Ehrlichia Antibody Rapid Test	Whole/Serum/Plasma
EHR-BG-ANA Ab	Ehrlichia/Babesia Gibsoni/Anaplasma Antibody Combo Rapid Test	Whole/Serum/Plasma
EHR-LYM-ANA-CHW	Ehrlichia/Lyme/Anaplasma/Henriworm Combo Rapid Test	Whole/Serum/Plasma
EHR-BG-ANA-CHW Ab	Ehrlichia/Babesia Gibsoni/Anaplasma/Heartworm Combo Rapid Test	Whole/Serum/Plasma
CPV-CCV-GIA-CRV Ag	CPV/CCV/GIA/CRV Combo Rapid Test	Feces
BC Ab		Whole/Serum/Plasma
	Babesia Canis Antibody Rapid Test	
CAV Ag	Canine Adenovirus Antigen Rapid Test	Oral, nasal, and eye secretor
CPIV Ag	Canine Parainfluenza Virus Ag Rapid Test	Oral, nasal, and eye secretion
HR-LYM-ANA-BAB Ab	Ehrlichia-Lyme-Anaplasma-Babesia Gibsoni Antibody Combo Rapid Test	Whole/Serum/Plasma
CDV-CAV-CIV-CPIV Ag	CPV/CAV/CIV/CPIV Antigen Combo Rapid Test	Oral, nasal, and eye secretor
CDV-CPV-CCV-GIA Ag	CDV/CPV/CCV/GIA Antigen Combo Rapid Test	Oral, nasal, and eye secretion
EHR-LSH Ab	Ehrlichia/Leishimania Antibody Combo Rapid Test	Whole/Serum/Plasma
HR-ANA-BG-BC Ab	EHR/ANA/BG/BC Antibody Combo Rapid Test	Whole/Serum/Plasma
Prog	Canine Progestrone Rapid Test	Serum/Plasma
HR+LSH+ANA+CHW	EHR Ah+LSH Ab+ANA Ab+CHW Ag Combo Rapid Test	Whole Blood/Serum/Plasma
CJBru Ag	Canine Brucella Antigen Rapid Test	Semen, vaginal secretions
CAV+CDV+CPV Ab	CAV+CDV+CPV Ab Combo Rapid Test	Whole Blood/Serum/Plasma
Angiostrongylus vasorum Ag	Angiostrongylus vasorum Ag Rapid Test	Whole Blood/Serum/Plasma
OB	Fecal Occult Blood Test	Feces
Anaplasma Ab+Dirofilaria Ag Ehrlichia Ab+Borrella Ab	Anaplasma Ab+Dirofilaria Ag+Ehrlichia Ab+Borrelia Ab Combo Rapid Test	Whole Blood/Serum/Plasma
Ovtomegalovirus Ag +Parvovirus Ag	Cylumegatownus Ag+Pervovirus Ag	
CPV-CCV-GIA-CRY Ag	Giandis+PARVO Vitrus-Cryposporidium+Corona virus Combo Rapid Test	Feces
EHR-ANA-Lyme-BAB-LSH	Ehlirichia+Anapiasma+Lyme disease+Babesia+Leishmania Combo Rapid Test	Whole Blood/Serun/Plasma



## Big Animal Rapid Test

#### 💮 Swine & Bovine & Sheep & Aviar

No.	Abbr.	Description	Uesing for	Sample
HFG1	TB MPB/T0.k 83	Sovine Tuberculosis Antibody (MPG70MO) Rapid test	Bovine	SerumiPlasma
HFGZ	PMDV-A	Anti-earthype A Food and Mouth Diesesse Antibody Repid Tost	Swine, Bovine, Sheep	SecuriPhoma
HF03	IBR Ab	Infectious Boxine Rhinotrachellis Anabody Ropid Test	Bovine	SerumPlasma
HF04	PAGS	Sovine Pregnancy Rapid Test	Bovine	Whole Blood /Serum/Plasma
HP95	5D (BCV-C pervers-ROTA-GRA-E coll K56)	Source 50 Arrigen Rapid Test (SCV-C, parvam/ROTA-GRA-E.coll K59)	Boyne	Peoes
HF06	BVDV Ab	Bovine Viral Diarrhee Virus Antibody Rapid Test	Bovine	Serum/Plosma
HF07	BCV Ag	Bovine Coronsvirus Antigen Rapid Test	Bovine	Feces
HFG8	ROTA Ag	Sovine Rolavirus Artigen Rapid Test	Bovine	Peoes
HF09	CE	Echinococcosis Antibody Rapid Test	Bovine Steep	Serum/Plosma
HF90	FMD NSP	Foot and Mouth Disease Virus Non-Structural Protein Antibody Rapid Test	Swine, Bovine, Sheep	SecuriPhona
HF11	BrutPS	Brucellows UPS antibody Rapid lest	Swine, Boyme, Sheep	Serum/Horma
HF12	PPR No	Pesta des Petits Ruminones Virus Arribady Rapid Tost	Sheep	SarumPloona
HF13	CSFV	Classical Swins Rever Virus Antigody Rapid Test	Swins	SecuriPlasma
HF14	FMDV-O	Anti-scrotype D Foot and Mouth Discusse Virus Antibody Rapid Test	Swine, Boyme, Sheep	SeumPlama
HF15	PRRSV	Potoine reproductive and respiratory syndrome virus Antibody Rapid Test	Swine	SerumiPlosma
HF95	PCV-2	Swine Porcine Circovirus 2 Genotypes Antibody Repid Test	Swine	SecuriPlanta
HF17	PRV-gB	Swine Pseudoniores Virus gB Antibody Rapid Test	Swine	SerumPlasma
HF18	TRI	Swine Trichinella antibody Rapid Test	Swine	SerumPlasma
HF19	PRV-gt	Swine Pseudorabies Visus glf. Antibody Rapid Test	Swine	SecuriPleana
HF20	ASFV Ab (P54)	African Swine Fever Virus Antibody (#54) Repid Test	Swine	SeumPlome
HF21	PEDV Ag	Poscine Epidenic Dianhoea Visus Antigen Rapid Test	Swine	Feces
HF22	TGEV Ag	Procine Transmissible Gastroenteritis Virus Antigen Rapid Test	Swine	Feces
HF23	NDV Ag	Newtonille Disease Virus Antigen Rapid Test	Avian	Traches, doses secretors, or leves
HF34	Salmonella Ag	Salmonella. Ag Rapid Test	Avian	Trachea, dozca secretions, or feces
HF25	ALV Ag	Asian Leukosis Wrus Ag Rapid Test	Avian	Traches, dosca secretions, or feces
HF28	AW Ag	Avian Influenza Vinus Ag Rispid Test	Avian	Traintee, doors secretors, or leses
HF27	BVDV Ag	Spvine hits#Diarrhea Virus Antigen Rapid Test	Bovine .	Ear tissue



Felina Plood Typing

## Blood Typing Test Kit - Feiline/Canine

"Fast, Accurate, Essential,"

## Feline Blood Typing Rapid Test

### Purpose:

Prevent acute hemolysis during transfusions.

### Key Facts:

Type B cats have strong anti-A antibodies; Type A cats have weak anti-B antibodies.

#### Risks:

Mismatched blood can cause vomiting, dyspnea, neurological signs, or death.

#### When to Test:

Pre-transfusion screening Anemia (pulse pressure <15%) Surgery (pulse pressure <20%) Hemolytic anemia Repeated transfusions

Sample: Whole blood (EDTA/heparin)

Time: 5-10 min





Result High accuracy, easy to read with naked eyes



















Canine Blood typin

Ringy bonis

## Canine DEA 1.1 Blood Typing Test

#### Detects:

DEA 1.1 status (positive/negative).

#### Use for:

Pre-transfusion screening Breeding risk assessment Preventing neonatal hemolysis

#### Result:

Rapid determination of DEA 1.1 phenotype

#### Importance:


Importance: DEA 1.1 is highly immunogenic; mismatched transfusion may cause:

Acute hemolytic reaction (RBC destruction within 12 hours)

Hemoglobinuria, hyperbilirubinemia


Delayed hemolytic reactions

Also triggered by: Pregnancy (≈25% of cases if mating with DEA 1.1+ male).













## **Bovine Pregnancy Rapid Test**

"Fast, Accurate, Essential,"

## Purpose:

Early detection of pregnancy in cattle (from 28 days post-breeding).

## Technology:

Double antibody sandwich immunoassay.

## Sample:

Serum, plasma, or whole blood.

## Advantages:

- 1. High accuracy, low cost
- 2. No device needed; visual results
- 3. Results in 20 minutes
- 4. After Bovine bleed or blood from Jugular vein or ear vein
- 5. Safe for embryos, reduces stress

## Application:

On-farm or lab use for all herd sizes.

## Result interpretation





## How to use?









<sup>\*</sup>Please interpret the result within 10 min. after the incubation.





## **Blood Collection Tube**

"Sharp. Sterile. Secure."









Volume: 0.2/0.25/0.5/1ml

Material: PP, chemical-resistant

Color-Coded Caps:

Green: Lithium heparin, chemistry Black: Sodium citrate, sedimentation

ellow: Gel & clot activator, serum immunology

Purple: EDTA, hematology (CBC) Red: No additive, general serum

Use: Veterinary micro blood collection for accurate lab diagnostics.



## Blood Collection Tube - Vacuum

"Sharp. Sterile. Secure."

## Available Sizes & Applications:

#### 2mL/3mL:

For small animals (cats, puppies, exotics)
— ideal for basic CBC or quick tests.

#### 5mL:

For medium animals or multi-test panels (e.g., biochemistry + immunity).

#### 10mL:

For large animals (horses, cattle) or blood bank storage.

## Features:

## Material Options:

Glass or PET and Centrifuge-safe

#### Vacuum Customization:

Altitude-adjusted vacuum levels (optional vacuum/non-vacuum)

## Versatility:

Suitable for a wide range of animal sizes and clinical needs









# **Blood pressure Cuff**

"Universal. Durable. Adaptable."

## Sizes:

1#-5# (3-15cm)

## Material:

Medical non-woven fabric

### Fits:

All animal sizes



#1: 3-6cm

Ej: turtle, weasel, mouse



#2: 4-8cm

Ej: cat, small dog, rabbit



Ej: dog, monkey



#4: 7-13cm

Ej: tiger, big dog

#5: 8-15cm

Ej: pig, lamp, ape





# Vet BP Monitor - SUPER Portable

"Accurate, Portable, Versatile,"









Compact & portable audible & visual

automatic shutdown

#### **PRODUCT FEATURES**

- . The measured values are accurate and sensitive.
- Optional SPO2 blood oxygen measurement.
- Standard with five blood pressure cuffs.
- Suitable for monitoring blood pressure, blood oxygen, and pulse rate of small, medium,
   and large animals such as cats and dogs, with simple operation and high cost-effectiveness.
- · Animal specific accessories, suitable for animals of different sizes.
- 4.3 inch high brightness color LCD display screen; Lightweight design for easy portability.
- High accuracy of blood oxygen saturation detection, supporting weak pulse detection function;
   Non invasive blood pressure measurement supports manual, automatic Continuous mode.

#### PRODUCT SPECIFICATION

NIBP		Systolic blood pressure	Diastolic blood pressure	Mean pressure		
	Large animals	40~270mmHg	10~215mmHg	20~230mmHg		
	Medium sized animal	40~200mmHg	10~150mmHg	20~165mmHg		
	Small animal	40~135mmHg	10~95mmHg	20~110mmHg		
	-	Static pressure range 0~300mmHg Static pressure accuracy±3mmHg				
	The second second					
Voltage		9-12VDC(with adpter)				
Max.in put power		70VA	/			
Fuse		F2AL 250V				
Battery capacity		2000mAh				
Temperature range		Working temperature: 5°C ~40°C				
		Transportation and storage temperature:-20°C ~50°C				
		Working humidity: 15%~80%				
Humidity range		Transportation and storage humidity: ≤93%				







# Vet BP Monitor - SUPER Portable

"Compact. Precise. Efficient."

## PRODUCT FEATURES

- Small size,easy to carry
- Slow inflation measurement. Equipped with pet-specific blood pressure cuff.
- Sensitive response time, easy communication, Rugged, easy to clean, safe and hygienic.
- Real time data storage, viewing historical data and trend graphs for review.
- · Precise measurement of blood pressure and pulse.
- Big data generates trend charts.



#### **PRODUCT FEATURES**

Model	V-E4D	NIBP	Pressure range 0mmHg~280mmHg	
Function	NIBP		Systolic blood pressure 60mmHg-255mmHg	
Power	Adapter(INPUTAC100-240V 50/60Hz		Diastolic blood pressure 30mmHg~195mmHg	
	OUTPUT DC5V 2A)	Working	Temperature:5°C~40°C	
Battery power supply	DC 3.7V	condition:	Humidity:15%~85%	
Product dimension	Approx.162*74*29mm		Atmospheric :80kPa~105kPa	
Product weight	Approx.250g	Transportation and	Temperature:-20°C~+55°C	
		storage condition:	Humidity:10%-95%	
			Atmospheric :50kPa~105kPa	



# Pet Wood's lamp

"Detect. Diagnose. Discover."

### Power:

15W, 4 LEDs

## Wavelength:

365nm (with black filter)

## Power:

1×18650 battery

## Design:

Aluminum alloy, compact (15.8×4.5×2.5cm), lightweight (150g)

### Use:

Detects fungal infections (e.g., ringworm), urine stains, and skin abnormalities in pets







18.6cm

7.5cm

## UV Wood's lamp - with Magnifier

"Portable, Precise, Safe,"

#### Use:

Handheld UV light for skin diagnosis in clinics

#### Features:

50mm lens, portable (182g), low RF interference

### Specs:

3.7V battery, works in 10°C-40°C





## Pet Thermometer - Advanced

"Accurate. Gentle. Reliable."



















Large LCD display, clear and easy to identify.





Two CR2032 batteries are provided for easy replacement.



10cm/3.9in

Product Name: Animal Body Temperat	ure Meter	After The End:Duzzer Prompt	
Measurement Time: About 20 Seconds   Display Meth		hod:LCD display(four and a haif digits)	
Power Supply:2*1.5V(AAA)	Resolution:±0.01°C(±0.01F)		
Operating Environment Temperature:	10°C~40°C(50	°F~104°F)	
Storage Temperature Of Instrument:0°C~50°C(32°F~122°F)			
Accuracy:±0.1°C 35°C~39°C(±2°F 95°F~102.2°F)			
Temperature Range:32.00°C~42.99°C(90.00°F~109.90°F)			

#### Metal Detector

Temperature measurement is quick, accurate, and easy to operate.







## Digital Thermometer - Soft Head Flexible Tip

"Accurate, Gentle, Reliable,"

## **Veterinary Thermometer**











123MM







Species	Normal temperature
Dog	100.2°F to 103.8°F (37.8°C to 39.9°C)
Cat	100.5°F to 102.5°F (38.0°C to 39.2°C)
Pig	101.6°F to 103.6°F (38.7°C to 39.8°C)
Cattle	100.4°F to 104.0°F (38.1°C to 39.9°C)
Sheep	100.9°F to 103.8°F (38.3°C to 39.9°C)
Horse	99.1°F to 100.8°F (37.3°C to 38.3°C)





## Digital Thermometer - Cover Of Flexible Tip

"Hygienic. Universal. Disposable."

100 Pack
Disposable Digital

#### Size:

98×26mm

## Packing:

100pcs/box

#### Features:

Disposable, hygienic, easy to apply

#### Use:

Universal fit for oral, rectal, and underarm thermometers

#### For:

Clinical and professional veterinary use only









# Long Gloves - Large Animal

"Soft. elastic. Environmental."





**EVA LONG GLOVE** SMOOTH/SUPER-SENSITIVE TYPE

It's made by 100% VIRGIN LDPE+EVA material.

COLOUR	SIZE	LENGTH	THICKNESS	WEIGHT	PACKING
	L	90 cm	22 mic	7.0 g	BOX/BAG
				olour vou want v	ve can supply.









# Long Gloves - Shoulder-Guard

"Durable, Sensitive, Protective,"



#### Use:

Bovine/equine rectal exams, AI, obstetrics

#### Material:

LDPE/EVA blend (strong, sensitive)

### Size:

85-93 cm length | Smooth surface

#### Colors:

Green, orange, red, blue

#### Pack:

50-100 pcs/box, bulk cartons available









# **Sharps Container**

"Safe, Secure, Durable,"

#### **Product Features:**

Eco-friendly PE material, anti-aging and high toughness Rotary seal cover prevents leakage Thickened design for durability









Rotate the Switch Clockwise




Unlocked, we put medical or biohazard waste into the container via the door.



after the switch shuted, the container CAN NOT be opened unless destroyed as so the sharps container is disposable.

## Opening Methods:

Cylindrical: Buckle lid to body (no rotation). Rotate red zone to lock permanently after use. Square: Flip red zone open or push cover forward. Close until click is heard after use.













## Sterilization Pouch - Self-Sealing

"Secure, Sterile, Smart."

#### Material:

Medical dialysis paper + PE/PET composite film








Easy hand-seal & peel, transparent view window Lead-free chemical indicator (color-changing arrows) Compatible with Steam/EO/Formaldehyde sterilization Triple-edge seal, anti-curling corner design

## Color change confirms sterilization:

brown (Steam), yellow (EO)











## Sterilization Reel - Heat-Sealing

"Versatile, Reliable, Sterile,"

#### Material:

Medical-grade paper (60/70gsm) + 5.2µm CPP/PET film (blue/green/white)



#### Sterilization Methods:

Compatible with Steam, EO, Gamma

#### Validity:

3 years


## Waterproof and tear resistant

Pre folding design and the adhesive strips


Can be folded easily and accurately.

#### Color Indicator:

Blue→Black (Steam), Pink→Yellow (EO)



200m



## Sterilizing reel roll



Product	Material	Size
		50mm*200m
	Made with high quality Medical grade paper and Composite membrance	55mm*200m
		75mm*200m
Sterilization Reel Roll		100mm*200m
		150mm*200m
		200mm*200m
		250mm*200m
		300mm*200m
		400mm*200m





## Sterilization Indicator Strips - Steam

"Reliable, Reactive, Recordable,"





- Widely used in hospitals, pharmaceuticals, food, healthproducts
- · Used for sterilization test
- Used as indicator whether the bag (or container) has beensterilized
- Sterilization indicator of medical equipments. The monitoring of the sterilization effect of the modelpressure steam sterilizer etc., whether it reaches the reference of sterilization.





#### Features:

Color change: Beige → Black (121°C/15min or 134°C/3min) Medical-grade paper with thermal-sensitive dye & adhesive Writeable, water-resistant, fade-resistant after sterilization



#### Sizes:

12/19/25mm × 50m (customizable)

## Package:

36/48/76 rolls/ctn (customizable)

## Usage:

For steam sterilization monitoring (pre-vacuum) Fixed packaging, process validation

## Storage:

Room temperature, dry, sealed, avoid light/acid/alkali



# Sterilization discoloration indication

Clear and accurate sterilization color change under steamsterilization.

# Medical gradepaper

Use medical profile paper asthe sub-stratecan be viscous, no residue, high adhesion, good flexibility simple. Cleaning and less pollution



# Peel without shatters

Unique design techniques, indicators will not be due tohigh pressure and self-shedding indication and peel withoutshatters





## Adjustable E-Collar

"Adjustable. Flexible. Protective."

# SIZE RECOMMENDATION

SIZE	Neck Girth	Depth	Recommended pets
1+	18.5"-20.8"	11"	large dogs, etc.
1#	16.1"-18.1"	8.6"	large dogs, etc.
2#	14.5"-16.1"	8"	Husky, Golden Retriever, Samoyed
2+	12.9"-14.5"	7.2"	Border Shepherd,etc.
3#	11.8"-13.3"	5.3"	Beijing Ba, Bichon, etc.
4#	10.6"-12.2"	5.3"	Poodle, Schnauzer, etc.
5#	8.6"-10.2"	5.3"	Pomeranian, Shih Tzu, etc.
6#	7"-8.2"	4.3"	Skinny cat etc
7#	6.2"-7.4"	3.5"	Kittens, Chihuahuas, etc.
8#	5.1"-6.2"	3.4"	kitten, puppies

Reminder: Due to the error of 1-3cm in manual measurement, please select the corresponding size specifications according to different sizes











## Pet E-Collar - Clear

"Adjustable. Flexible. Protective."



#### Material:

100% polypropylene (flexible & transparent)

### Fit:

Easy to apply, 9 sizes for all dogs

#### Use:

Post-surgery, wound care, skin diseases, hotspots

### Benefits:

Clear view, comfortable recovery





	epth Neck girth
Neck girth	Depth

		The second second	The state of the s		
Size	Neck girth(cm)	Depth(cm)	Open size(cm)		
1#	74	40.5	130*65		
2#	68	35.5	114*58.5		
3#	58	30	97.5*50		
4#	50	25	81.5*42		
5#	38	20	67.5*34.5		
6#	38	15	55*30		
7#	33	12.5	45*24		
8#	25	10	37*20		
9#	25	7.5	32*17.5		



## Pet E-Collar - Soft Fabric

"Soft. Safe. Customizable."

Sizes:

S, M

Material:

Non-woven fabric

Colors:


Blue, Dark blue (custom colors available)

For:

Dogs & cats | Prevents biting and licking of wounds













## Velux Super Collar

"Comfortable, Protective, Durable."

The VETLUX SUPER COLLAR redefines pet recovery wear.

Designed for superior comfort and protection, it allows for natural movement while effectively preventing pets from licking or biting their wounds.

Its durable, yet lightweight construction ensures a stress-free recovery period.





Available in precise sizes to ensure a perfect fit for all breeds:







## **Pet Lint Roller**

"Effective. Convenient. Essential."

#### Versatile Use:

Perfect for clothing, upholstery, blankets, and car seats.

## Compact & Portable:

Lightweight design ideal for home, travel, or on-the-go touch-ups.

## **Eco-Friendly Options:**

Refill packs available to reduce waste (60 or 90 sheets).



Easily lifts and traps pet hair, lint, and dust in seconds.

## Tear-Off Sheets:

Disposable sheets for mess-free cleaning. Available in 60 or 90 sheets per roll.

## Options:

Starter Kit: 1 Roller + 1 Refill (60 or 90 sheets)
Refill Pack: 60 or 90 sheets (compatible with standard lint rollers)











## Pet Leash

"Strong. Light. Reliable."

## Dimensions:

5 ft (152.4 cm) length × 0.4 inch (1.0 cm) diameter

## Material:

Abrasion-resistant, colorfast polyethylene (PE) braid

## Strength:

Supports over 300 lbs tensile force

#### Features:

Welded metal ring (2.0 cm), lightweight (60 g), easy to clean

#### Colors:

Blue/White, Blue/Black, Green/Black, Red/Black

#### Use:

Ideal for veterinary and professional handling, tangle-resistant, high visibility







Good Quality Hook

Other

# Pet Leash - with Zinc Hook

"Durable, Secure, Stylish."

## Reinforced Zinc-Alloy Hook:

Durable, bite-resistant, and rust-proof for maximum security.

## Reinforced Stitching:

Machine-stitched edging on handle for exceptional strength and longevity.

## Multiple Color Options:

Available in various colors to match your style.







## Pet First Aid kit

"Ready. Comprehensive. Essential."

## 101 PCS

Name: cross border pet first aid kit

Material: EVA+Oxford fabric

Size: 26\*19\*8cm

Weight: 0.7kg

Color: Red

Product category: 26 types 101 pcs set

Application: Emergency wound

treatment for cats and dogs

Features: Convenient to carry, Various

accessories









Large Capacity

Portable Handle

Customized Logo







INSIDE













## **Acupuncture Needles**

"Sterile, Versatile, Precise,"

## Material:

304 Stainless Steel with Copper Handle

## Available Sizes (Diameter × Length):

0.20 × 13mm 0.20 × 25mm

0.22 × 25mm 0.22 × 40mm

0.25 × 25mm 0.25 × 40mm 0.25 × 50mm 0.25 × 60mm

0.30 × 25mm 0.30 × 40mm 0.30 × 50mm 0.30 × 75mm

 $0.35 \times 75$ mm







# Dental scale tip

"Precise, Durable, Professional,"

## Model:

G1/G2/G3/G4/G5/P1

#### Material:

Precision-crafted from medical-grade stainless steel for exceptional durability and corrosion resistance.

## Application:

Designed for professional dental cleaning, polishing, and filling procedures.

## Compatibility:


Ideal for use with standard air polisher units in clinical settings.

## Origin:

Manufactured in HeBei, China, with strict quality control.

## Features:

High-performance design ensures efficient plaque removal and smooth polishing.













## **Dental Bur Stand**

"Precise, Durable, Professional,"

Imported 304 Stainless Steel: High corrosion resistance and durability

Textured Comfort Grip: Secure non-slip handling, even in wet conditions

Heat & Chemical Resistant: Suitable for repeated sterilization (autoclave/chemical)

Solid One-Piece Construction: No loose parts, enhanced stability and longevity

Hole Options: 3, 4, 5, 6, 8, or 10 holes available to suit various surgical needs











## Handpiece Sterilization Tray - Low-Speed

"Precise. Durable. Professional."

## Versatile Storage:

Holds low-speed burs, polishing brushes, prophylaxis angles, and scaler tips securely.



## Autoclavable:

Withstands high-pressure steam sterilization up to 135°C.

## Flexible Organization:

Available in 4, 8, 12, or 18-hole configurations.

#### **Durable & Efficient:**

Designed for repeated clinical use and optimal infection control.







## **Endoscopic Linear Cutter**





Model	Shaft Length	Articulation	Silicone Handle
VLELCS01	60mm	Yes	No
VLELCS02	160mm	Yes	No
VLELCS03	260mm	Yes	No
VLELCS04	60mm	Yes	Yes
VLELCS05	160mm	Yes	Yes
VLELCS06	260mm	Yes	Yes







## **Endoscopic Reloader**

"Reload, Reposition, Resume."





There are various specifications of components available tomeet different needs based on the height and effectivecutting length of the matching nail.







Precise and efficient reduction of surgical trauma, helpingpatients recover faster.







- · 11 adjusting degrees
- 45° maximum left and right
- Convenient and accurate access to target tissue and blood vessels
- · Easy access to narrow spaces
- Reduce damage and risk of tissue separation
- · More convenient to hold tissue
- Improve power of tissue closure
- Ensure the firing effect of each pin is consistent





# Endoscopic Reloader

"Reload. Reposition. Resume."

Model	Description	Color	STAPLE SIZE (INNER TO OUTER ROW)	Stapler length
VLELCS10	Curved Tip Articulating Vascular Reload	Gray	2-2-2mm	30mm
VLELCS11	Curved Tip Articulating Vascular/ Medium Reload with Unequal height staplers Technology	Tan	2-2.5-3mm	30mm
VLELCS12	Curved Tip Articulating Voscular/ Medium Reload with Unequal height staplers Technology	Tan	2-2.5-3mm	45mm
VLELCS13	Curved Tip Articulating Medium/Thick Reload with Unequal height staplers Technology	Purple	3-3.5-4mm	30mm
VLELCS14	Curved Tip Articulating Medium/Thick Reload with Unequal height staplers Technology	Purple	3-3.5-4mm	45mm
VLELCS15	Curved Tip Articulating Medium/Thick Reload with Unequal height staplers Technology	Purple	3-3.5-4mm	60mm
VLELCS16	Curved Tip Articulating Extra Thick Reload with Unequal height staplers Technology	Black	4-4.5-5mm	45mm
VLELCS17	Curved Tip Articulating Extra Thick Reload with Unequal height staplers Technology	Black	4-4.5-5mm	60mm





# **Endoscopic Trocar**

"Access. Seal. Introduce."



Reference	Shaft Diameter(mm)	Shaft Length (mm)
VLELCS07	5	100
VLELCS08	10	100
VLELCS09	12	100







# BOOST YOUR CURE



Vetmed - Headquarter Address: Saurine Baeva 1, Almaty, Kazakhstan Tel: (7) 777 217 00 99 E-mail: Max@vetmed.com

Vetlux Vet Consumable Limited - Subsidiary UNIT 1307, BEVERLEY COMM CENTRE, 87-105 CHATHAM RD,SOUTH, TSIM SHA TSUI,HONG KONG Tel: (86) 13145820826 / (86) 17673469281 E-mail: Jack@vetluxgroup.com / Kenny@vetluxgrou.com